


# SIDEKICK: The Mobile Bucky

A Mobile Positioning Companion for Mobile DR


With **powered vertical and horizontal travel**, plus **powered tilt and pivot motions**, let the Sidekick do all the work for you. No more awkward adjustments being made within the personal space of your patient, or valuable time being spent fine tuning without visibility of your frame of reference. Position your Sidekick with the press of a button!




- Wireless remote controlled motion
- Actuator-driven vertical and horizontal motions
- Actuator-driven panel tilt and pivot motions
- Coarse and fine adjustment settings
- Ergonomic panel loading and unloading
- Sanitizer and disinfectant resistant surface
- Locking casters



## Footprint





## Telescoping Ranges of Motion




VERTICAL  
+  
HORIZONTAL  
TELESCOPING

## Panel Head Ranges of Motion



## Omni Wheel Add-On



TechnoAide