

1. System Overview

1.1. Application

- Abdomen
- Gynecology
- Obstetrics
- Cardiology
- Small parts
- MSK
- Urology
- Vascular
- Pediatric

1.2. Transducer Types

- Convex array
- Linear array
- Endocavity curved array
- Endocavity linear array

1.3. Imaging Modes

- B-mode
- M-mode
- Pulsed Wave Doppler

1.4. Imaging Technique & Function

- Harmonic Imaging
- Multi-Beam-Forming
- Speckle Reduction Imaging (eSRI)
- Synthetic Receiving Aperture
- Pan Zoom
- Needle Guide

1.5. Display Modes

- B
- B+B
- 4B
- M
- B+M
- B+PW

1.6. System Language Support

- Software : English, Chinese, German, French, Italian, Russian, Spanish, Polish, Portuguese

- Keyboard input: English, Chinese, Russian
- Control panel overlay: English, Chinese
- User Manual: English, Chinese

1.7. Standard Features

- USB Disk (only for international marketing)
- Probe cable holder
- Dustproof cloth

1.8. Optional Features

- Transducers
- Printers
- Battery (only for international marketing)
- Hard disk
- Needle Guide Bracket Kit
- Mobile trolley
- Footswitch
- Hand carried bag
- Luxury Hand carried bag
- DICOM 3.0

2. Physical Specification

2.1. System Architecture

- Physical Channels: 24
- System dynamic range: 0-158dB
- Beam forming: Dual beam
- Memory: 504MB
- Hard drive: 320GB(optional)
- Operation System: ucos

2.2. Dimension and weight

- Height: 32cm
- Width: 22cm
- Depth: 33cm
- Weight: 7.1kg

2.3. Monitor

- 12.1" TFT-LCD monitor
- Resolution: 1024 x 768
- Imaging field size: 640*512
- Video out size: 1024 * 680
- Capture size: BMP, JPG, AVI, DCM 1024 * 768; FRM, CIN 640 * 512

- View angle: Up 80°, Down 80°, Left 80°, Right 80°
- Brightness, Contrast and Color Temp adjustable
- Built-in stereo speaker

2.4. Transducer port and holder

- 2 active ports
- 2 holders
- 1 Coupling holder

2.5. Electrical Power

- Voltage: 110V-240VAC
- Frequency: 50/60 Hz
- Power: 74.4w

2.6. Battery

- Rechargeable lithium ion battery
- Capacity: 5000mAh
- Removable
- Approximately 60 minutes of typical ultrasound exam use
- Max charging time: 8 hours

2.7. Environmental operating requirements

- Ambient temperature: 5° to 40°C
- Relative Humidity: 25%~80% (no condensation)
- Atmospheric pressure: 86kPa-106kPa

2.8. Environmental storage requirements

- Ambient temperature: -20° to 55°C
- Relative Humidity: 25%~93% (no condensation)
- Atmospheric pressure: 70kPa-106kPa

3. User Interface

3.1. Control Panel

- Interactive back-lighting
- Hard Keys provides tactile feedback
- Physical trackball
- 8 segment TGC sliders
- Physical keyboard

3.2. System boot-up

- Boot up from complete shut-down in about 30 sec
- Shut-down in about 3 sec
- Recovery from screen saver in about 3 sec

3.3. Comments

- Arrow
- Block move and delete for separate blocks of text
- Support physical keyboard for text input
- 150 comments in comments library
- 15 User customizable comments per a library

3.4. BodyMark

- Up to 130 Body Mark graphics in library

3.5. Screen Information

- EDAN logo
- Hospital name
- Date
- Time
- Patient ID
- Patient Name
- Transducer model
- Preset name
- Mechanical index (MI)
- Thermal Index (TI)
- Imaging parameters
- Gray Scale bar
- Depth Scale

***Not all the items are listed in here, please refer to the User Manual.**

4. Imaging Parameters

4.1. B-Mode

- Zoom: 7 levels, x1.44, x1.96, x2.56, x4.0, x5.76, x9.0, x16.0 (available on live state); 3 levels: x1.78, x4.0, x16.0 (available on freeze state)
- Depth: 1.9- 32.4cm

- Frequency: 2.0-15.0MHz(3 fundamental & 2 harmonic frequencies)
- eSRI: 0-3
- Rejection: 0-7
- Scan Angle: 0-3
- Gain: 0-130dB, 2dB/step
- TGC: 8 segments
- Dynamic range: 30-150 dB, 4dB/step
- Scan Density: H/M/L
- Max Frame rate: 270 f/s
- Map: 0-14
- Pseudo Color: 6 types
- Frame Persist: 0-7
- Focus position: 0-15
- Focus Number: 1-4
- GAC: 0-7
- IP: 0-7
- SRA: 0-1
- H Reverse: On/Off
- V Reverse: On/Off
- 90° rotation : 0/90/180/270
- B/W Invert: On/Off
- Display format: single(B), dual(B+B), Quad(4B): B, 2B, 4B

4.2. M-Mode

- Sweep speed: 0-3
- Line Average: 0-7
- Gray Map: 0-14
- Pseudo Color: 6 types
- Gain: 0-130dB, 2dB/step
- Frequency: 2.0-15.0MHz(3 fundamental frequencies)
- Dynamic range: 30-150 dB, step 4dB/step
- B/M Display: L/R, Full

4.3. Pulsed Wave Doppler

- Frequency: 2.5-6.5MHz, 2 levels
- PRF: 13 levels
- Gain: 0-130dB, 2dB/step
- Dynamic range: 30-150dB, 4dB/step
- Wall filter: 0-3
- Sweep speed: 0-3
- Baseline: 0-6
- Correct Angle: 15° -165°, 1°/step

- Steer: 0°, ±10°(available on linear transducers)
- Invert: Up/Down
- Volume: 0-7
- Pseudo color: 6 types
- Sample Volume: 16 levels, 0.5-20 mm
- D Rejection: 0-7

5. Cine Review and Post-Processing

5.1. Cine Review

- Frame by frame manual review
- Auto playback
- Start frame and end frame are selectable for cine loop review
- Maximum cine memory is up to:
256 frames for B mode
23s for M mode
10s for PW Doppler mode

5.2. Post-Processing

- B Mode: zoom, pseudo color, Gray map
- M Mode: pseudo color, Gray map
- PW: pseudo color, Gray map

***Not available the stored images and clips in**

Review

- FRM and CINE file support measurement, comments, bodymark

6. Imaging Storage and Exam

Database

6.1. Imaging Storage

- 504MB for data storage 504MB
- 320GB hard drive
- Storage up to approximately 224 static images as BMP format(memory); storage up to approximately 145635 static images as BMP format (hard drive)
- Maximum clip is up to:
256 frames for B mode

23s for M mode
10s for PW Doppler mode

6.2. File Management

- Support exam storage temporarily without patient information
- Support image files query
- Support delete, rename image files
- Support review image files of current exam or prior exam
- Support store images as BMP, JPG, FRM, AVI, CIN or DCM format
- Support export images to a USB disk

7. Connectivity

- DICOM Storage:
Verify SCP
Static image store SCU DCM
Removable media
Manual-transfer on demand
- 2 USB Ports
- Video out:
VGA
Video: PAL/NTSC
S-video: PAL/NTSC
- Footswitch port
- Remote port
- Ethernet

8. Preset

- Application:
Abdomen
Obstetric
Gynecology
Small Parts
Urology
Vascular
Cardiac
Pediatric
- Presets
Abdomen
Abd Difficult

Aorta
Obstetric
Gynecology
Endovaginal
Small Parts
Urology
Vascular
Cardiac
Pediatric

9. Peripheral & Accession

9.1. Printer

- Black/white Digital/Analog Video printer
SONY UP-X898MD
MITSUBISHI P93W_Z
- Color Analog video printer
SONY UP-25MD
SONY UP-D25MD
- Graph/text printer
HP LaserJet Pro 400 M401d402D 403D
HP LaserJet 1510
HP Deskjet 1010
HP DeskJet Ink Advantage Ultra 2029
HP DeskJet 1112

9.2. Needle Guide Bracket

- **BGK-CR60**
 - Focus Depth:44mm
 - Angle:45°
- **BGK-LA43**
 - Focus Depth:26.5mm
 - Angle:43°
- **BGK-CR10UA**
 - Focus Depth: 250mm
 - Angle: 3°
- **BGK-LA70**
 - Focus Depth:42.5mm
 - Angle:44°
- **BGK-MCR10**
 - Focus Depth:20mm
 - Angle:35°
- **BGK-EL40**

- Focus Depth:23mm
- Angle:0°

10. Measurement and Report

10.1. General Measurement

- **B Mode**
 - Distance
 - Cir/Area
 - Volume: 2-Axis, 3-Axis, 3-Axis(LWH)
 - Ratio
 - Angle
 - %Stenosis:Distance, Area
 - Histogram
- **M Mode**
 - Distance
 - Time
 - Slope
 - Heart Rate
- **Doppler**
 - Velocity
 - Heart Rate
 - Time
 - Acceleration
 - RI
 - Auto: PS, ED, TAMAX, RI, PI, S/D
 - Trace Direction: Above, below, Dual
 - Trace Sensitivity+
 - Trace Sensitivity-

10.2. Application Measurement

- **Gynecology**
 - B mode**
 - UT:Length, Width, Height, UT
 - Endo
 - OV-Vol(L/R): Length, Width, Height, OV-Vol
 - FO(L/R, Number:4): Length, Width
 - CX-L
 - UT-L/CX-L
 - PW mode**
 - Velocity, L UT A, R UT A, L OV A, R OV A, Trace Direction, Trace Sensitivity+, Trace Sensitivity-

- **OB**
 - B mode-OB-1**
 - GS , CRL, BPD, HC, AC, FL, AFI, EFW
 - FBP
 - Growth Curve
 - EDC
 - B mode-OB-2**
 - TAD, APAD, CER, FTA, HUM, OFD, THD, NT
 - FBP
 - EDC
 - PW mode**
 - Velocity, FHR, Umb A, MCA, Fetal AO, Desc.AO, Placent A, Ductus V, Trace Direction, Trace Sensitivity+, Trace Sensitivity-
 - M mode**
 - FHR, Time, Slope
- **Cardiac**
 - M mode**
 - LV: TEICHHOLZ(LVIDd, LVIDs, ET, HR, EDV, EDS, SV, CO, EF, FS, SI, CI, MVCF, BSA), CUBE(LVIDd, LVIDs, ET, HR, EDV, EDS, SV, CO, EF, FS, SI, CI, MVCF, BSA)
 - Mitral Valve: EF Slope, ACV, A/E, Valve Volume
 - Aorta: LAD/AOD, Valve Volume
 - Heart Rate
 - LVET: LVET, AVSA
 - LVMW: LVPWd, IVSTd, LVIDd, LVMW
 - B mode**
 - RV
 - LV: S-P Ellipse (LVLd, LVALd, LVLs, LVALs, EDV, ESV, SV, CO, EF, SI, CI, BSA), B-P Ellipse (LVALd, LVAMd, LVIDd, LVALs, LVAMs, LVIDs, EDV, ESV, SV, CO, EF, SI, CI, BSA), Bullet (LVAMd, LVLd, LVAMs, LVLs, EDV, ESV, SV, CO, EF, SI, CI, BSA), Mod. Simpson(LVAMd, LVLd, LVAPd, LVAMs, LVLs, LVAPs, EDV, ESV, SV, CO, EF, SI, CI, BSA)
 - RV
 - PA
- **Small Part**
 - B mode**

- L.THY-V: Length, Width, Height, Volume
- R.THY-V: Length, Width, Height, Volume
- Isthmus

● **Urology**

B mode

- BLV: Length, Width, Height, Volume
- RUV: Length, Width, Height, Volume
- Prostate Vol: Length, Width, Height, Volume, PPSA, PSAD

● **Pediatric**

B mode

- HIP(α , β)

● **Vascular**

PW mode

- Velocity, CCA, ICA, ECA, Vert A, Upper, Lower, Trace Direction, Trace Sensitivity+, Trace Sensitivity-

***For more measurement information, please refer to the User Manual.**

10.3. Report

- Worksheet
- Diagnostic
- Configure whether print image in report
- Export as PDF format

11. Transducers

C361-2

- Imaging Format: convex array
- Number of Elements: 80
- Convex Radius: 60 mm
- FOV: 60°
- Bandwidth: 2.0-6.0MHz
- Fundamental Frequency: 2.5 MHz, 3.5MHz, 4.5MHz
- Harmonic Frequency: H5.0MHz, H5.4MHz
- Doppler Frequency: 2.5MHz, 3.0MHz
- Depth: 19-324mm
- Frame Rate(18cm, Full of FOV): max 103 f/s
- PW velocity: max 2.01m/s($\pm 60^\circ$)
- Applications: Abdomen, OB, Gynecology, Urology
- Needle Guide Bracket: BGK-CR60(14G, 18G,

20G, 22G)

L743-2

- Imaging Format: general linear array
- Number of Elements: 96
- Footprint: 40 mm
- Bandwidth: 5.0-10.0MHz
- Fundamental Frequency: 6.5 MHz, 7.5MHz, 8.5MHz
- Harmonic Frequency: 9.0 MHz, 9.4MHz
- Doppler Frequency: 5.5MHz, 6.5MHz
- Depth: 19-117mm
- PW velocity : max 1.14m/s($\pm 60^\circ$)
- Applications: SMP, Vascular
- Needle Guide Bracket: BGK-LA43(14G, 18G, 20G, 22G)

L761-2

- Imaging Format: general linear array
- Number of Elements: 80
- Footprint: 60mm
- Bandwidth: 5.0-10.0MHz
- Fundamental Frequency: 6.5 MHz, 7.5MHz, 8.5MHz
- Harmonic Frequency: 9.0 MHz, 9.4MHz
- Doppler Frequency: 5.5MHz, 6.5MHz
- Depth: 29-117mm
- PW velocity : max 1.14m/s($\pm 60^\circ$)
- Applications: SMP, Vascular
- Needle Guide Bracket: BGK-LA70(14G, 18G, 20G, 22G)

E611-2

- Imaging Format: endocavity micro convex array
- Number of elements: 80
- Convex Radius: 10 mm
- FOV: 155°
- Bandwidth: 5.0-8.0MHz
- Fundamental Frequency: 5.5 MHz, 6.5MHz, 7.5MHz
- Harmonic Frequency: 9.0 MHz, 9.4MHz
- Doppler Frequency: 5.0MHz, 6.0MHz
- Depth: 19-157mm
- Frame Rate(10cm, Full of FOV): max168 f/s
- PW velocity: max 0.56m/s($\pm 60^\circ$)
- Applications: OB, Gynecology, Urology

- Needle Guide Bracket: BGK-CR10UA(16G)

C611-2

- Imaging Format: micro convex array
- Number of elements: 80
- Convex Radius: 10 mm
- FOV: 128°
- Bandwidth: 5.0-8.0MHz
- Fundamental Frequency: 5.5 MHz, 6.5MHz, 7.5MHz
- Harmonic Frequency: 9.0 MHz, 9.4MHz
- Doppler Frequency: 5.0MHz, 6.0MHz
- Depth: 29-157mm
- Frame Rate(10cm, Full of FOV): max 168 f/s
- PW velocity : max 0.79m/s($\pm 60^\circ$)
- Applications: Pediatric, Pediatric Cardiac
- Needle Guide Bracket: BGK-MCR10(14G, 18G, 20G, 22G)

E741-2

- Imaging Format: endocavitylinear array
- Number of elements: 80
- Footprint: 40 mm
- Bandwidth: 5.0-10.0MHz
- Fundamental Frequency: 6.5 MHz, 7.5MHz, 8.5MHz
- Harmonic Frequency: 9.0 MHz, 9.4MHz
- Doppler Frequency: 5.5MHz, 6.5MHz
- Depth: 19-117mm
- PW velocity : max 0.69m/s($\pm 60^\circ$)
- Applications: Urology
- Needle Guide Bracket: BGK-EL40(16G, 18G)

Engineering

- EN ISO 14971: Medical Device Risk Management
- ISO 10993: Medical Device Biocompatibility
- NEMA UD 2: Output Measurement for Diagnostic Ultrasound Equipment

12. Regulatory approvals

- FDA Class II Device
- CE/MDD Class IIa
- IEC 60601-1: Medical Equipment Safety
- IEC 60601-1-2: Medical Device Electromagnetic Safety
- IEC 60601-2-37: Ultrasonic Medical Equipment Safety
- IEC 62133: Battery Safety
- IEC 62304: Medical Device Software Life-cycle Process
- IEC 62366: Medical Device Usability