

# OptiWell™ Combo Pack

*Deep Well Plates + Tip Combs for KingFisher™*



- Our most popular plate and comb in one pack
- An ideal pair for automated research platforms
- Made from medical-grade polypropylene
- Compatible with KingFisher™ platforms - Apex, Duo Prime, Flex and Presto

OptiWell™ Tip and Comb sterile combo packs are a convenient package built for compatibility with KingFisher™ and Accuris IsoPure™ automation platforms.



The included 96-well plate features square wells with conical v-bottoms, a working well volume of 2.2mL and printed alphanumeric.



OptiWell™ tip combs are made from a low-binding, medical grade polypropylene and are ideal for applications in DNA extraction / quantification, viral RNA extraction, immunoprecipitation and more.

OptiWell™ Sterile Combo Packs are packaged in cases of 50 coupled pairs of Tip Combs and Deep Well Plates.

| Cat. No.                                       | Description  | Qty |
|--|--|-----|
| D3096-KFC                                      | OptiWell™ Tip Comb & Plate Combo Pack - Includes:<br>96-Well Plate w/ Square Wells, PP, Conical V-Bottom, 2.2mL,<br>96-Tip Comb, PP, Conical V-Bottom, Sterile, 10 Bags of 5 pairs | 50  |
| <i>Also available for individual purchase:</i> |  |     |
| D3096-39S                                      | Deep Well Plate, 96 Square Well, 2.2mL, V-Bottom, 44.2mm Height, Sterile   | 50  |
| D3340-96S                                      | Tip Comb, 96-Tip, V-Bottom, PP, for KingFisher™, Sterile   | 100 |

\*KingFisher™ is a trademark of Thermo Fisher Scientific