

Classic Express, Flat Top, Treatment Table with Full Shelf

MODEL 200-E

- New, ultra-durable, PVC-laminated, solid, hardwood legs with levelers
- Easy assembly with pre-threaded inserts and hanger bolts in the legs
- Full shelf in easy-clean laminate
- Seamless rounded corner top
- Paper dispenser included
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Bright Maple finish
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 400 lbs. load capacity under normal use (181 kg)



Maple and Black

SPECIFICATIONS

200-E27	72" L x 27" W x 31" H
200-E30	72" L x 30" W x 31" H



Maple and Royal Blue

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.