

Classic Express, Flat Top, Treatment Table with Drawers

MODEL 1300-E

- New, ultra-durable, PVC-laminated, solid, hardwood legs with levelers
- H-brace construction
- Laminate shelf with two easy clean, plastic drawers
- Legs have pre-inserted carriage bolts for fast easy assembly
- Seamless, rounded corner top
- 2" of firm, foam padding
- Paper dispenser included
- Available with 27" and 30" tops
- Bright Maple finish
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 400 lbs. load capacity under normal use (181 kg)



Maple and Black



Maple and Royal Blue

SPECIFICATIONS

1300-E27	72" L x 27" W x 31" H
1300-E30	72" L x 30" W x 31" H



California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2