

Classic Express, Treatment Table with Drawers

MODEL 1013-E

- New, ultra-durable, PVC-laminated, solid, hardwood legs with levelers
- Easy assembly with pre-threaded inserts and hanger bolts in the legs
- Laminate shelf, two easy clean, plastic drawers and adjustable backrest
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Drawers can be positioned at either end, and from either side
- Fast and easy assembly
- Seamless rounded corner top
- Paper dispenser included
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Bright Maple finish
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 400 lbs. load capacity under normal use (181 kg)

SPECIFICATIONS

1013-E27 72" L x 27" W x 31" H

1013-E30 72" L x 30" W x 31" H



Maple and Black



Maple and Royal Blue

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1