

Classic Express, Treatment Table with Shelving

MODEL 1030-E

- New, ultra-durable, PVC-laminated, solid, hardwood legs with levelers
- Easy assembly with pre-threaded inserts and hanger bolts in the legs
- Laminate, tiered shelf, storage compartments and adjustable backrest
- Adjustable backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Fast and easy assembly
- Seamless rounded corner top
- Paper dispenser included
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Bright Maple finish
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 400 lbs. load capacity under normal use (181 kg)

SPECIFICATIONS

1030-E27	72" L x 27" W x 31" H
1030-E30	72" L x 30" W x 31" H



California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2



Maple and Black



Maple and Royal Blue