

Alpha Express, Treatment Table with Drawers

MODEL 3013-E

- Gray, powder-coated, all-steel frame, gray laminate shelf, two, easy-clean, plastic drawers with gray drawer fronts and an adjustable backrest
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 450 lbs. load capacity under normal use (204 kg)



Black



Royal Blue

SPECIFICATIONS

3013-E27	72" L x 27" W x 31" H
3013-E30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1