

Alpha Express, Treatment Table with Shelving

MODEL 3030-E

- Gray, powder-coated, all-steel frame, gray laminate shelving and an adjustable backrest
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 450 lbs. load capacity under normal use (204 kg)

SPECIFICATIONS

3030-E27	72" L x 27" W x 31" H
3030-E30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2



Black



Royal Blue