


Pure•Amp™ PCR Plates

Clear and White 384-Well PCR Plates

- 🔍 Ultra-thin walled tube wells
- 🔍 Alphanumeric reference grids
- 🔍 RNase, DNase, DNA, and PCR inhibitor free
- 🔍 Compatible with most thermal cyclers:
 - Bio-Rad™
 - ABI™
 - Roche®
 - Etc.



Introducing the newest additions to the Pure•Amp™ family of PCR plates!



Pure•Amp™ 384-well PCR plates are fully skirted with the appropriate notches to work with most automation platforms. They are also offered in natural/clear and opaque white for greater fluorescence reflectivity.

Pure•Amp™ 384-Well PCR Plate – Platform Compatibility Chart

P9384-C	P9384-W
Applied Biosystems (ABI™) – 384 Well Variants:	
ViiA7 / 7900HT / Quantstudio 5, 6, 7, and 12K / 3130xL, 3500xL, 3730xL	
Bio-Rad™ – 384 Well Variants:	
CFX Opus 384	CFX Opus 384
Roche® – 384 Well Variants:	
	LightCycler 480I & 480II (384 Well)
Analytik Jena – 384 Well Variants:	
	qTower ³ 84

ABI™ and its respective brand names are trademarks of Thermo Fisher Scientific | Bio-Rad™ and its respective brand names are trademarks of Bio-Rad Laboratories, Inc. Roche® and its respective brand names are trademarks of Roche | Analytik Jena and its respective brand names are a trademark of Analytik Jena, an Endress+Hauser Company






P9384-C - Natural Color
ABI™ + Bio-Rad™ Type Plate



P9384-W - White Color
Roche + Bio-Rad™ Type Plate



Cat. No.	Description	Qty
Pure•Amp™ 384 Well PCR Plates		
P9384-C	PCR Plates, 384 well, 30µL, full skirt, natural, ABI® type, one notch - A24	50
P9384-W	PCR Plates, 384 well, 30µL, full skirt, white, Roche® type, two notch - A24/P24	50
Pure•Amp™ 96×0.2mL Well PCR Plates		
P9602-N-T	Sep-R-8™ Breakable PCR Plates, Standard Well, 96×0.2mL, Non-Skirted, Breaks into Strips of 8	50
P9602-N	PCR Plates, Standard 96×0.2mL Well, Non-Skirted	50
P9602-NW	PCR Plates, Standard 96×0.2mL well, Non-Skirted, White	50
P9602-S	PCR Plates, 96×0.2mL Well, Semi Skirted, Natural	50
P9602-SRR	PCR Plates, 96×0.2mL Well, Semi Skirted w/ Raised Rim, Natural (ABI Type)	50
Pure•Amp™ 96×0.1mL Well PCR Plates		
P9601-N	PCR Plates, 96×0.1mL Well (Low Profile/Fast), Non-Skirted	50
P9601-NW	PCR Plates, 96×0.1mL (Low Profile/Fast), Non-Skirted	50
P9601-S	PCR Plates, 96×0.1mL (Low Profile/Fast), Semi Skirted	50
P9601-SRR	PCR Plates, 96×0.1mL Raised Rim (Low Profile/Fast), Semi Skirted	50
P9601-SW	PCR Plates, 96×0.1mL (Low Profile/Fast), Semi Skirted, White	50