

U-Leg Couch



MODEL 3600

- Dual frame construction with separate base & headrest sections
- Non-adjustable wedge headrest
- Chrome-plated, steel legs
- Paper dispenser included
- 2" firm foam padding
- Heavyweight, knitbacked upholstery
- 250 lbs. load capacity under normal use

SPECIFICATIONS


3600-24	72" L x 24" W x 18" H
3600-27	72" L x 27" W x 18" H

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Options



040 Paper Holder



016 Adjustable Wedge Headrest

Upholstery Colors

