

# CONVERTIBLE STAIRCASES

## STANDARD FEATURES:

- Load capacity: 350 lbs. (Maple Laminate) & 450 lbs. (Maple & Birch Wood)
- Stairs can be configured as a straight or corner unit
- Platform: 30"W x 30"D x 24"H
- (4) 30"W x 9"D x 6"H steps and (6) 30"W x 9"D x 4"H steps
- Handrails are made from solid hardwood (shipped ready to assemble)
- Uprights are made from 2" steel with a black powder coat finish
- Stairs are available in wood or laminate finishes
- Anti-slip safety tread on platform and all steps for increased safety
- 9' Ceiling height required



PTST-005-WHR

## STAIRCASE CONFIGURATION:

### PTST-005-WHR

Platform: 30"W x 30"D x 24"H

Stairs: 102"L x 30"W x 24"H

9"D Steps

Overall including Handrails (straight): 102"L x 34"W x 59"H

Overall including Handrails (corner): 78"L x 60"W x 59"H

### PTST-005-WHR-ADA

Platform: 36"W x 36"D x 24"H

Stairs: 125"L x 36"W x 24"H

11"D Steps

Overall including Handrails (straight): 125"L x 40"W x 59"H

Overall including Handrails (corner): 94"L x 72"W x 59"H

## STANDARD WOOD CHOICES:



Birch



Maple

## STANDARD LAMINATE CHOICES:



Maple