

# DHH-10

DIGITAL HOME HEALTH SCALE



**RICE LAKE**<sup>®</sup>  
WEIGHING SYSTEMS





# DHH-10

## DIGITAL HOME HEALTH SCALE

Rice Lake's DHH-10 digital home health scale is designed to efficiently weigh up to 440 pounds (200 kilograms). Users can alternate between pound and kilogram weighing with a simple press of a button. Featuring a bright LCD display and clear weighing surface, the DHH-10 is battery operated and indicates to users when power is low. An optional DHH-10BT Bluetooth® model can be synced with the AccuroFit application on Apple® and Android® mobile devices for personalized health monitoring.

### STANDARD FEATURES

---

- Dimensions: 15.4 in x 11.8 in x 1 in (W x L x H)
- Capacity: 440 lb x 0.2 lb (200 kg x 0.1 kg)
- 8 mm tempered safety glass weighing surface
- Battery operation, (2) AAA alkaline batteries included
- Low battery indication
- Display lb or kg
- Large LCD display

### WARRANTY

---

- One-year limited warranty

