

QuikSip™

QuikSip™ bottle-top aspirator

The QuikSip™ is a convenient bottle-top aspirator that removes supernatants, cell culture media, and other biological liquids from dishes, flasks, culture bottles, centrifugation tubes, multiwell plates, and more. It is small, portable, and doesn't require an external vacuum or power source.

- **Moves Easily Between Labs:** The QuikSip™ works without any external vacuum source, allowing it to be easily transferred between workstations. Each plunger stroke powers up to 25mL of aspiration.
- **Provides Fingertip Aspiration Control:** The QuikSip™ includes a unique thumb-operated controller that is compatible with disposable pipette tips, micropipettes, and glass Pasteur pipettes.
- **Adapts to Sterile Environments:** All wetted parts are autoclavable (except piston/cylinder cartridge). Sterile piston/cylinders are available for use in sterile applications.
- **Fits Most Reagent Containers:** Mounts on most common reagent bottles with its 45mm thread and 45/33mm and 45/38mm adapters (included). Other thread adapters available.

The Economy Aspirator for Biological Fluids.



Technical Data

Operating limits and exclusions

Liquids unsuitable for use with the QuikSip™ bottle-top aspirator include, but are not limited to:

- concentrated or oxidizing acids
- non-polar solvents, such as hydrocarbons and halogenated hydrocarbons
- hydrofluoric acid (HF)
- those that attack PP, PE, EPDM and silicone
- temperature of instrument and reagents must be between 15° and 40°C (59° and 104°F)
- piston/cylinder cartridge and pump housing are not autoclavable

Items supplied

QuikSip™ aspirator includes:

- QuikSip™ controller with 200cm (78in.) of silicone tubing
- Spare piston/cylinder cartridge
- Two PP thread adapters 45/33mm and 45/38mm. See page 40 for other adapters
- Operating manual
- One-year warranty



QuikSip™ aspirator shown with optional autoclavable 8-channel manifold.

ORDERING INFORMATION

Description	Cat. No.
QuikSip™ bottle-top aspirator	4724150
8-channel manifold; PP, autoclavable, pack of 1	704526
25mL replacement piston/cylinder cartridge, pack of 3	704504

Bottles & Adapters

See detailed listing on page 40-41.