



The Performalift has always offered the ultimate in comfort with the ability to tilt the table up at the head end to relieve sinus pressure or tilt it in the opposite direction with a salon or backrest top to create a zero gravity chair position.

Now available in ADA compliant height or high height range.

PERFORMA LIFT SERIES



Flat Top



Lift Assist Backrest



Electric Salon



Zero Gravity

The unique twin tower lift system creates a Zero Gravity chair capability insuring the ultimate in comfort. The tilting function is helpful for many treatments and offers greater comfort for people who experience sinus pressure when lying flat.



Programmable Hand Control

Now standard, this control has two user programmable memory positions and the ability to change two settings at once with a single button. The simple interface makes positioning clients fast and the new electronics are quiet and smooth in operation.



Integrated Heated Top (Optional)

Our new system is imbedded in the foam, creating a more even heat, and has three settings. Designed originally to be used in car seats and hospital beds, the heating element automatically shuts off after 10 hours providing peace of mind for safety and energy conservation.

SPECIFICATIONS

- The unique twin tower lift system
- Tilttable top for greater comfort
- Programmable memory positions
- Replaceable AeroCel™ padding
- Now available in ADA compliant height range and high height

TABLE SPECIFICATIONS

Lift System	Electric, Twin Tower
Controls	Foot Control + Hand Control with memory positions
Lifting Capacity	550 lbs. (249.5 kg)
Fabric	TerraTouch - Medical Grade, ISO 10993
Finish	Water Based Stain
Padding	4" (10.2 cm) or 6" (15.2 cm) AeroCel™ (Optional)
Height Range 4" (10.2 cm) pad	20" - 34" (51 cm - 86 cm) or 19" - 30" (48 cm - 76 cm)
Tilt	+/- 18°
Width	29" (73.6 cm), 31" (78.7 cm), 33" (83.8 cm)
Length	73" (185.4 cm)
Options (must be ordered with table)	SoViVibroacoustic System, Integrated Heated Top, ABC System, Casters
Warranty	2 years parts, 1 years labor
Safety Listings	cETLus Listed product 3034177 to UL 962, CE Marked

Accessories	Adjustable Side Arm Rests, Boiance Shelf Support, Boiance Face Rest Cushion, Quick Lock Face Rest
-------------	---

WOOD FINISHES



TERRATOUGH COLORS


Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing. Please refer to the product cleaning instructions for a list of approved cleaners and disinfectants.




OPTIONS



Patented ABC System (Adjustable Body Comfort) is an incredible improvement in comfort for women.



SoVi vibroacoustic system is designed to deliver all the latest forms of healing modalities via a very flexible integrated delivery system.



The Fuzion heating system is embedded in the foam, creating a more even heat.


ACCESSORIES



The Boiance Facerest Pad with its patented water sphere design offers unmatched comfort when coupled with our highly adjustable QuickLock Facerest.



The Boiance Shelf Support creates the perfect place for clients arms. The shelf can also be used alone as a head support.



The adjustable side armrests are versatile and comfortable and feature a wide range of positions. (Electric Salon only)