

# DETECTO®

# 7550

## CONNECTED WALL-MOUNT FOLD-UP WHEELCHAIR SCALE

 Bluetooth™

 Wi-Fi  DS smart



Stylish gray color and unique fold-up storage allow the 7550 scale to integrate into any healthcare environment.

# WHEELCHAIR SCALE BLENDS INTO ANY HEALTHCARE ENVIRONMENT



Install the 7550 Right in High-Traffic Areas and Fold Up When Done Weighing



- Retains patient dignity by weighing in open commons areas instead of backrooms or closets
- The fold-up non-weighing mode allows you to store the scale in high-traffic areas and work zones without using up valuable floor space
- The wall-mount scale stores easily out of the way and aesthetically blends into any environment for a more efficient workflow
- Slim 7-inch / 18-cm profile takes up minimal space when installed in a medical facility hallway, clinical office, or patient room
- 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity
- Bluetooth 5.0 / Wi-Fi included standard for wireless connectivity with EMR/EHR systems
- Sturdy steel construction for long-lasting durability over a lifetime of use

Integrated full-color, touchscreen LCD with wheelchair tare weight mode to remove wheelchair weights and factor patient weight



**DETECTO**



Slim  
7.3-in/  
18.5-cm  
Compact  
Profile

# 7550

## FOLD-UP, WALL-MOUNT DIGITAL WHEELCHAIR SCALE

### Mounts Securely to the Wall for Storage

- 36 in W x 33 in D / 91 cm W x 84 cm D Weighing Platform Provides Ample Room
- Wheelchair Tare Entry Can be Done Before or After the Patient is on the Scale
- Accommodates Some Sizes of Power Wheelchairs
- May be Used on Carpet or Hard Floor Surfaces
- Wi-Fi and Bluetooth 5.0 Wireless Connectivity
- 6 C Batteries (not included) or AC Adapter (included on some models)
- Body Mass Index Calculation
- Units Locking into LB or KG
- Two-Way Easy Access Ramps
- Auto Weight Lock Feature
- Gas Spring Platform Lifting Mechanism
- StableSENSE Adjustable Filtering
- Textured Non-Skid Mat on Ramp and Platform
- Compact Vertical Fold-Up Platform Conserves Space
- Pushbutton and Keypad Tare for Wheelchair Weights

### 8-Digit Numeric Patient ID

DETECTO's 7550 incorporates an iPhone-like touchscreen LCD that allows you to enter up to 8-digit patient IDs that can be sent along with patient weight data to EMR systems via Wi-Fi or Bluetooth 5.0.



Detecto's model 7550 wheelchair scale features a unique fold-up storage design that saves space when not in use. The 7550 easily accommodates most wheelchairs, and yet requires very little floor space when folded up. The slim 7-inch / 18-cm profile takes up minimal space when installed in a medical facility hallway, clinical office, long-term care center, or nursing home.

When you're ready to weigh, simply lower the platform and the gas spring lifting mechanism will gently lower it into place. The 7550's 1.5-inch/38-mm high platform is very low to the ground for easy patient access. This makes the process simple and minimal effort required in order to roll a patient on and off the weighing platform.





# 7550

## CONNECTED WALL-MOUNT FOLD-UP WHEELCHAIR SCALE



	7550	7550-AC	7550KGEU
<b>Capacity</b>	1,000 lb x 0.2 lb / 450 kg x 0.1 kg		300 kg x 0.1 kg
<b>AC Adapter</b>			●
<b>OIML EC Class III</b>			M
<b>Platform Size (Flat Portion)</b>	36 in W x 33 in D / 91 cm W x 84 cm D		
<b>Width</b>	50.3 in / 127.6 cm		
<b>Depth</b>	39.3 in / 99.7 cm (with platform extended); 7.3 in / 18.4 cm (without platform extended)		
<b>Overall Height</b>	44.0 in / 111.7 cm		
<b>Display Type</b>	3.2 in (diagonal) full color TFT display, 320 x 240 resolution		
<b>Number of Characters</b>	Weight: 5 digits, 0.8 in (20 mm) high Height: 4 digits, 0.35 in (9 mm) high BMI: 3 digits, 0.35 in (9 mm) high		
<b>Keys</b>	Mechanical switch type, Power, Zero, Lock/Release, ID, Cycle Display Mode, Menu, Print		
<b>Power Requirements</b>	Six C size Alkaline, Ni-Cad, or NiMH batteries (not included) Domestic models: Optional 100 to 240V AC 50/60Hz 12V DC 1A wall plug-in UL/CSA listed AC power adapter (Cardinal part number 6800-1045) International models: Optional 100 to 240V AC 50/60Hz 12V DC 1A Multi-Pin Input AC power adapter (Cardinal part number 6800-1047)		
<b>Operating Environment</b>	Operating Temperature Range: 14 to 104 °F (-10 to +40 °C), Humidity: 0 to 90% non-condensing		
<b>Communication Interfaces</b>	Wi-Fi and Bluetooth 5.0		
<b>Construction</b>	Steel		
<b>Weighing Mechanism</b>	Electronic load cells		
<b>Lifting Mechanism</b>	Gas springs		
<b>Net Weight</b>	180 lb / 82 kg	181 lb / 83 kg	180 lb / 82 kg
<b>Shipping Weight</b>	238 lb / 108 kg	239 lb / 109 kg	238 lb / 108 kg

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.



**SOLD BY:**



[www.universalmedicalinc.com](http://www.universalmedicalinc.com)