

# SO-LOW

ENVIRONMENTAL EQUIPMENT COMPANY

SETTING TEMPERATURE STANDARDS IN EVERY LAB

EST. 1959



## ULTRA LOW FREEZER & LAB REFRIGERATION EQUIPMENT PRODUCT CATALOG

MEDICAL RESEARCH | INDUSTRIAL USE | COLD STORAGE | LAB TESTING | LIFE SCIENCES



MEDICAL RESEARCH | PHARMACEUTICAL | COLD STORAGE | CLINICAL TRIALS | LIFE SCIENCES

### Advanced Touchscreen Interface

The Platinum Series Ultra Low Freezers come with an advanced PID control system. Freezers are controlled using a 4.3" full color touch screen display. Control is eye-level door mounted on upright style models. Temperature recording, downloadable via USB drive (CSV format), embedded SCADA, and LAN features.


### Electrical Approval

So-low Freezers are (ETL) listed and conform to UL Standard 61010-1.

So-low Freezers are Canada/CSA Standard C 22.2 No. 51010-1.

### Nationwide Service Network

Technical problems are serviced by a Nation-Wide-Network of service companies.



MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	DEFROST TYPE	UNIT WEIGHT
--------------	------------	--------	-------------------	------------------------------------	------------------------------------	--------------	-------------

#### UPRIGHT VERTICAL TO -85°C

NU85-13	13	370	-40°C to -85°C	20" x 22" x 51.25"	36" x 33.5" x 79.5"	MANUAL	675 LBS
NU85-18	18	505	-40°C to -85°C	23" x 26" x 51.25"	36.5" x 36.5" x 79.5"	MANUAL	725 LBS
NU85-22	22	626	-40°C to -85°C	28.5" x 26" x 51.25"	42" x 36.5" x 79.5"	MANUAL	800 LBS
NU85-25	25	710	-40°C to -85°C	34" x 25" x 51.25"	49.5" x 36.5" x 79.5"	MANUAL	900 LBS
NU80-28	28	792	-40°C to -80°C	38" x 25" x 51.25"	53" x 36.5" x 79.5"	MANUAL	950 LBS
NU80-30	31	843	-40°C to -80°C	42" x 25" x 51.25"	59.5" x 36.5" x 79.5"	MANUAL	975 LBS

#### CHEST HORIZONTAL TO -85°C

NC85-2	2	56	-40°C to -85°C	17" x 14" x 14"	25" x 26" x 38.5"	MANUAL	300 LBS
NC85-3	3	83	-40°C to -85°C	24" x 12" x 18"	34.5" x 26" x 47"	MANUAL	350 LBS
NC85-5	5	142	-40°C to -85°C	30" x 18" x 16"	40.5" x 32" x 47"	MANUAL	425 LBS
NC85-9	9	255	-40°C to -85°C	47" x 16" x 20"	58" x 29.5" x 47"	MANUAL	500 LBS
NC85-12	12	340	-40°C to -85°C	55" x 19" x 20"	66" x 32.5" x 48"	MANUAL	650 LBS
NC85-14	14	400	-40°C to -85°C	62" x 20" x 20"	73" x 34" x 48"	MANUAL	725 LBS
NC85-17	17	480	-40°C to -85°C	69" x 20" x 21"	81" x 34.5" x 49"	MANUAL	775 LBS
NC80-21	21	594	-40°C to -80°C	81" x 21.5" x 21"	91" x 36" x 49"	MANUAL	825 LBS
NC85-22	22	626	-40°C to -85°C	67" x 19" x 30"	98" x 36" x 44"	MANUAL	850 LBS
NC80-27	27	766	-40°C to -80°C	67" x 23" x 31"	98" x 36" x 44"	MANUAL	900 LBS