


Casco C-Matt™ OR Table Pads

Casco C-Matt™ OR Table Pads are designed with both the patient and the surgeon in mind by providing a surface that keeps the patient safe, clean, and comfortable while operating.



Casco C-Matt™ OR Table Pads come in a variety of sizes to best accommodate both your facility's and patients' needs, coming in thicknesses of 2, 3, and 4 inches. OR Table Pads can be purchased as a complete set, or as individual component parts. Our OR Table Pads are manufactured to fit popular Skytron and Amsco-Steris OR Table Models. Armboard Pads are universal fit.

Casco C-Matt™ OR Table Pads are made from a breathable, bi-directional stretch polycarbonate fabric with a high density, antimicrobial, hypoallergenic medical grade foam core. Our OR Table Pads are also durable, easy to clean, stain resistant, and RF-sealed to keep fluid from penetrating the product. There is also a non-skid material on the bottom to keep the product and patient secure.

