



Casco C-Matt™ Positioning Wedges

Casco C-Matt™ Positioning Wedges provide therapeutic solutions to a wide variety of users, ranging from infants to the elderly, through a variety of products designed to fit the patient's needs.



Casco C-Matt™ Positioning Wedges come in a variety of makes, models, shapes, and sizes to best accommodate both your facility's and patients' needs, coming in angles of 15, 20, 35, 45, and 55 degrees, all of which are made to support the head, hips, knees, spine, cervix, and more.



Casco C-Matt™ Positioning Wedges are made from a breathable, bi-directional stretch polyurethane fabric with a high density, antimicrobial, hypoallergenic medical grade foam core. Our Positioning Wedges are also durable, easy to clean, stain resistant, and RF-sealed to keep fluid from penetrating the product.



Product Specifications

Product Name	Part # GTIN UPC
15 Degree Positioning Wedge - 9.5" Depth x 10.5" Wide x 2.55" Height	H12619 00815073029243 815073029243
20 Degree Positioning Wedge - 8.5" Depth x 11" Wide x 3.09" Height	H12620 00815073029250 815073029250
35/55 Degree Positioning Wedge - 9.5" Depth x 14" Wide x 6.65" Height	H12621 00815073029267 815073029267
45 Degree Positioning Wedge - 7" Depth x 24" Wide x 7" Height	H12622 00815073029274 815073029274
Tall Knee Bolster Positioning Wedge - 17.675" Depth x 21" Wide x 10" Height	H12623 00815073029281 815073029281
Knee Bolster Positioning Wedge -12.5" Depth x 21" Wide x 7.25" Height	H12642 00815073029458 815073029458
Between the Knee Support Positioning Wedge - 7.9" Depth x 10.2" Wide x 6.1" Height	H12643 00815073029465 815073029465
Head Support Donut Positioning Wedge -9" diameter	H12644 00815073029472 815073029472
Head Cradle Positioning Wedge -8" Depth x 9" Wide x 4.25" Height	H12645 00815073029489 815073029489