

C-Matt Safety Securement Seating System

Extraordinary support and comfort for people in wheelchairs.

Key Benefits

Casco Prevention Wheelchair Securement Pads provide an additional means of enhancing comfort in a wheelchair by providing a breathable, easy to clean, adjustable, and durable means of positioning the user in an upright and comfortable manner..



Product Features

Zoned Medical Grade Foam

- Density formulated for comfort and durability which works well in a long term period of sitting,

Polycarbonate surface fabric

- Absorbs body heat and perspiration, dispersing it away from the body, keeping the patient cool and dry reducing the incidence of skin break down.

The Casco Prevention Wheelchair Securement Pad is uniquely designed to redistribute pressure to prevent wound development and to promote healing of existing ones


C-Matt Safety Securement Seating System

Identification Codes

CASCO Part Number	H12762
GTIN	00810125100765
UPC	810125100765



- The securement buckles are rated at a 600 lbs. break strength.
- The webbing has an elongation 2% stretch rating and a break strength well over 1000lbs.
- The pad is easy to install and remove and can be cleaned using a mild dish detergent or antibacterial cleaner.
- The harness is easy to adjust to different sizes.



BASED ON A 2 HOUR TIME TEST
(MMHG-millimeters of mercury)

Footnote: This product does not cause entrapment – meaning that the patient is not confined to the product, and can leave at any time thanks to the five point harness.