



POLTEX

CLEANING AND DISINFECTING INSTRUCTIONS FOR PETG PRODUCTS

Cleaning:

Wash PETG products with a mild soap or detergent and plenty of lukewarm water. Use a soft cloth, applying only light pressure. Rinse with clear water and dry by blotting with a damp cloth or chamois.

Easy to clean with any typical hospital cleaning solution, even 10% bleach.

Disinfecting:

Use a soft cloth with any of the below solutions:

10% bleach solution

70% Isopropyl Alcohol

Up to 96% Ethyl Alcohol

Up to 28% Hydrogen Peroxide

Clorox or other disinfecting wipes

If your disinfectant is not listed above and have a question of reactivity with PETG, test an inconspicuous area before cleaning the whole product.