

**WIDE UNIVERSAL ARM BOARD**

**PRODUCT INFORMATION**

**COPYRIGHT ©2025. ALL RIGHTS RESERVED**



**PHILLIPS SAFETY**  
PRODUCTS INC.

## PRODUCT INFORMATION

WIDE UNIVERSAL ARM BOARD | AB-11.5x31.5-WIDE



The Wide Universal Arm Board is designed to provide both patient comfort and clinician efficiency during procedures that require reliable arm positioning. Its expanded width accommodates both the patient's arm and essential procedural supplies, giving staff a stable, convenient workspace without the need for additional accessories.

Crafted from durable HDPE, this arm board is radiolucent and MR-safe, ensuring clear imaging compatibility across surgical, radiology, and interventional environments. Its tapered design allows the board to slide easily beneath the patient's shoulder, supporting smoother insertion and promoting effective venous distention—especially valuable during IV placement, injections, or contrast administration.

The smooth HDPE surface cleans easily, supporting strict hygiene and infection control protocols.

The Phillips Safety Wide Universal Arm Board is an excellent choice for facilities that need a versatile, easy-to-position support surface that improves workflow and enhances patient stability throughout a wide range of procedures.

## SPECIFICATION

---

**COLOR:** White

**MATERIAL:** HDPE

**SIZE:** 31.5" x 11.5" x 3/8"

**FEATURES:** Tapered design;

Radiolucent; MRI Safe.

**PRODUCT INFORMATION**  
 WIDE UNIVERSAL ARM BOARD | AB-11.5x31.5-WIDE

