

# PHENOLIC UNIVERSAL ARM BOARD

PRODUCT INFORMATION

COPYRIGHT ©2025. ALL RIGHTS RESERVED



## PRODUCT INFORMATION

PHENOLIC UNIVERSAL ARM BOARD | AB-36x6



The Phenolic Universal Arm Board is designed to provide stable, reliable support during procedures requiring the injection of contrast materials or sustained arm positioning. Made from high-quality phenolic, this arm board delivers exceptional durability while remaining radiolucent and MR-safe, ensuring clear imaging without interference across a wide range of modalities.

Its lollipop-style design makes it ideal for surgical tables without side rails, offering quick placement and secure positioning where traditional rail-mounted boards aren't compatible. The smooth phenolic surface is easy to clean and disinfect, supporting stringent infection control protocols in surgical, radiology, and interventional environments.

With a generous size of 36" x 6", the board provides excellent support for a wide variety of patient arm positions, helping clinicians maintain proper access and ergonomics during procedures.

The Phillips Safety Phenolic Universal Arm Board is a dependable, versatile solution for operating rooms, cath labs, and imaging suites seeking a long-lasting, easy-to-use support surface for patient care.

## SPECIFICATION

---

**COLOR:** Black


**MATERIAL:** Phenolic

**SIZE:** 36" x 6" x 1/4"

**FEATURES:** Lollipop-style design;

Radiolucent; MRI Safe.

**PRODUCT INFORMATION**  
 PHENOLIC UNIVERSAL ARM BOARD | AB-36x6



PRODUCT INFORMATION

PROJECT

Universal 6 x 36 Phenolic Board

CREATED FOR: PHILLIPS-SAFETY

APPROVED:	SIZE	WEIGHT	DWG NO	REV
CHECKED:	A3		Universal 6 x 36 Phenolic Board	
CREATED: Jerry Verden / etx3d.com		12/5/2025	SCALE 1:5	