


MODEL NU80-30

30 CUBIC FOOT ULTRA LOW TEMPERATURE FREEZER


| | WIDTH | DEPTH | HEIGHT | WIDTH | DEPTH | HEIGHT |
|---------------------|---------|-------|----------|--------|--------|--------|
| INTERIOR DIMENSIONS | 42 IN | 25 IN | 51.25 IN | 106 CM | 63 CM | 130 CM |
| EXTERIOR DIMENSIONS | 60.5 IN | 37 IN | 79.5 IN | 153 CM | 93 CM | 201 CM |
| CRATED DIMENSIONS | 68 IN | 45 IN | 86 IN | 172 CM | 114 CM | 218 CM |

NOTE THAT SHELVES ARE FIXED AND CANNOT BE REMOVED OR ADJUSTED UNLESS SPECIAL ORDERED AS ADJUSTABLE.

SIDE VIEW



FRONT VIEW



*TO PASS THROUGH DOORWAYS, THE DOOR CAN OPEN A FULL 180°, TO REDUCE THE UNIT EXTERIOR DEPTH TO 33"

SO-LOW



NET WEIGHT : 975 LBS.
SHIPPING WEIGHT: 1325 LBS.

SHIPPING OPTIONS

- DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)
- DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)
- INSIDE DELIVERY (WHITE GLOVE DELIVERY)