



C-Matt™ Comfort

A Basic Prevention Surface to Meet Your Facility's Health, Safety and Patient Comfort Needs



Entry Level Prevention Performance

Multi-zoned Foam Core for Patient Comfort and Safety

Leak Proof Top Cover

Head to Foot Reversible

Optional Fire Barrier for Regulatory Compliance

Patient Weight Limit: 350lb

Top Quality Replacement Surfaces from Casco

Surface Foam Core

- *High density, anti-microbial, medical grade foam*
- *Medium-firm support with convoluted topper for added comfort*
- *Firm foam perimeter for patient safety*

Regulatory Requirements

- *Meets CAL117*
- *Optional fire barrier meets CAL129/FED1633*

Replaceable Cover

- *Anti-microbial and hypo-allergenic*
- *Easy to clean, stain resistant and durable*
- *Zipper-flap design to protect foam core*
- *Clamshell design for easy maintenance*

Standard (6 x 36 x 80 / 84) and Custom Sizes Available

5 year Warranty

Clamshell cover design keeps the top and bottom cover connected



Convoluted foam topper and firm perimeter



Quality replacement surfaces at a 30-50% cost savings.

