



# C-Matt™ Comfort Ultima

*A High Performing Prevention Surface to Meet Your Facility's Skin Safety and Comfort Needs*



## The Ultimate in Patient Comfort

### Conforming Viscoelastic Memory Foam Topper

*Works with your patient's unique body shape to redistribute pressure*

### Breathable, Bi-Directional Stretch Top Cover

### Shear Liner for Patient Comfort

### Optional Fire Barrier for Regulatory Compliance

**Patient Weight Limit: 350lb**

## Top Quality Replacement Surfaces from Casco

### Surface Foam Core

- *Viscoelastic memory foam topper conforms to patient's body for optimal comfort*
- *High density, anti-microbial, medical grade foam*
- *Firm perimeter and base for patient safety*

### Regulatory Requirements


- *Meets CAL117*
- *Optional fire barrier meets CAL129/FED1633*

### Replaceable Cover


- *Anti-microbial and hypo-allergenic*
- *Easy to clean, stain resistant and durable*
- *Patented zipper-flap design to protect foam core*
- *Clamshell design for easy maintenance*

**Standard (6 x 36 x 80 / 84) and Custom Sizes Available**

Clamshell cover design keeps the top and bottom cover connected



Specialized foam zones for optimal pressure redistribution



Quality replacement surfaces at a 30-50% cost savings

