

SO-LOW

ENVIRONMENTAL EQUIPMENT COMPANY

SETTING TEMPERATURE STANDARDS IN EVERY LAB

EST. 1959



ULTRA LOW FREEZER & LAB REFRIGERATION EQUIPMENT PRODUCT CATALOG

MEDICAL RESEARCH | INDUSTRIAL USE | COLD STORAGE | LAB TESTING | LIFE SCIENCES



MEDICAL RESEARCH | PHARMACEUTICAL | COLD STORAGE | CLINICAL TRIALS | LIFE SCIENCES

So-Low's Superior Construction

Manufactured at our facility in Cincinnati Ohio USA, So-Low's Low Temperature Freezers to -40°C are produced with high quality craftsmanship, heavy-duty materials, and over 60 years of experience. These freezers use a hermetic motor compressor refrigeration system, providing durable and reliable environmental conditions for your product's vital temperature requirements.

Electrical Approval

So-low Freezers are (ETL) listed and conform to UL Standard 61010-1.
Canada/CSA Standard C 22.2 No. 51010-1.

Nationwide Service Network

Technical problems are serviced by a Nation-Wide-Network of service companies.



MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	DEFROST TYPE	UNIT WEIGHT
--------------	------------	--------	-------------------	------------------------------------	------------------------------------	--------------	-------------

UPRIGHT VERTICAL TO -40°C

U40-13	13	370	0°C to -40°C	20" x 22" x 51.25"	36" x 33.5" x 79.5"	MANUAL	675 LBS
U40-18	18	505	0°C to -40°C	23" x 26" x 51.25"	36.5" x 36.5" x 79.5"	MANUAL	725 LBS
U40-22	22	626	0°C to -40°C	28.5" x 26" x 51.25"	42" x 36.5" x 79.5"	MANUAL	800 LBS
U40-25	25	710	0°C to -40°C	34" x 25" x 51.25"	49.5" x 36.5" x 79.5"	MANUAL	900 LBS
U40-28	28	792	0°C to -40°C	38" x 25" x 51.25"	53" x 36.5" x 79.5"	MANUAL	950 LBS

CHEST HORIZONTAL TO -40°C

C40-2	2	56	0°C to -40°C	17" x 14" x 14"	25" x 26" x 38.5"	MANUAL	300 LBS
C40-3	3	83	0°C to -40°C	24" x 12" x 18"	34.5" x 26" x 47"	MANUAL	350 LBS
C40-5	5	42	0°C to -40°C	30" x 18" x 16"	40.5" x 32" x 47"	MANUAL	425 LBS
C40-9	9	255	0°C to -40°C	47" x 16" x 20"	58" x 29.5" x 47"	MANUAL	500 LBS
C40-12	12	340	0°C to -40°C	55" x 19" x 20"	66" x 32.5" x 48"	MANUAL	650 LBS
C40-14	14	400	0°C to -40°C	62" x 20" x 20"	73" x 34" x 48"	MANUAL	725 LBS
C40-17	17	480	0°C to -40°C	69" x 20" x 21"	81" x 34.5" x 49"	MANUAL	775 LBS
C40-21	21	594	0°C to -40°C	81" x 21.5" x 21"	91" x 36" x 49"	MANUAL	825 LBS
C40-22	22	626	0°C to -40°C	67" x 19" x 30"	98" x 36" x 44"	MANUAL	850 LBS
C40-27	27	766	0°C to -40°C	67" x 23" x 31"	98" x 36" x 44"	MANUAL	900 LBS