

# MODEL C40-27


## 27 CUBIC FOOT ULTRA LOW TEMPERATURE FREEZER

|                     | WIDTH  | DEPTH | HEIGHT   | WIDTH  | DEPTH  | HEIGHT |
|---------------------|--------|-------|----------|--------|--------|--------|
| INTERIOR DIMENSIONS | 67 IN  | 23 IN | 31.25 IN | 170 CM | 57 CM  | 79 CM  |
| EXTERIOR DIMENSIONS | 98 IN  | 36 IN | 43.5 IN  | 248 CM | 91 CM  | 110 CM |
| CRATED DIMENSIONS   | 101 IN | 43 IN | 53 IN    | 256 CM | 109 CM | 134 CM |

NOTE THAT THE FREEZER LID CAN OPEN UP TO APPROXIMATELY 75° FROM HORIZONTAL POSITION

### SIDE VIEW

### FRONT VIEW



# SO-LOW



**NET WEIGHT : 900 LBS.**  
**SHIPPING WEIGHT: 1125 LBS.**

### SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY ( REQUIRES FORKLIFT OR PALLET JACK )  
 INSIDE / WHITE GLOVE DELIVERY ( REQUIRES A SEMI-TRUCK RECEIVING DOCK )