

## MODEL C40-17X

### EXPLOSION PROOF CHEST FREEZER 0°C TO -40°C

<b>Capacity [ Cu. Ft. ]</b>	<b>17 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>480 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>0°C to -40°C</b>
<b>Temperature Range [ °F ]</b>	<b>32°F to -40°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	81.25" x 37.5" x 48.5"
Interior [ Inches ]	69" x 20" x 21"
Exterior [ Centimeters ]	206 x 95 x 123
Interior [ Centimeters ]	175 x 50 x 53
Refrigeration System	Single Stage System
Compressors	1, Hermetic Compressor
Compressor Container	Explosion Proof
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	HCFC & CFC Free
Exterior	Powder Coated Steel
Interior	Zinc Coated Galvanized Steel
Insulation	Polyurethane (Ecomate) Foam
Additional Insulation	Beaded Foam Sub-Lids
Lockable Lid	Lockable with Padlock
Casters	8, Swivel with Locks
Storage	Optional Inventory Racks
Shipping Weight	975 Lbs.

#### Electrical

Standard Voltage	115 Volts, 60 Hertz, 1 Phase
Amps / Run Amps / Plug ( Optional Voltages )	20A   7.3A   NEMA 5-20 208 Volts, 230 Volts
Power Configuration	Hardwired
Certification(s)	UL Rated
Certification File Number	SA10711
Certification Classification	Class I, Division I & II, Groups C & D

#### Controller

Display Type	Explosion Proof, Digital
--------------	--------------------------

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 1 Year Parts & Labor