

PRODUCT INFORMATION

MOBILE LEADED BARRIER WITH 36" X 24" WINDOW



The Mobile Leaded Barrier with 36" x 24" Window is a portable, full-body radiation barrier. It provides a barrier of 2.0mm lead equivalence steel, and features a large 36" W x 24" H viewing panel with 2.0mm lead equivalency glass and 0.5mm lead equivalence acrylic. In addition, the Phillips Safety Mobile Leaded Barrier has steel legs with heavy duty casters and a locking caster system for ease of mobility. With an overall size of 39.5" W x 75" H, it is recommended that two people move the barrier together at all times. Plus, it weights 220 lbs.

SPECIFICATION

WINDOW SIZE: 24" x 36"

BARRIER SIZE: 39.5" W x 75" H

BARRIER MATERIAL: Metal

PROTECTION LEVEL: Lead Glass 2.0 PbEq,
Lead Acrylic 0.5mm PbEq

FEATURES: Viewing Panel and Locking caster system

PRODUCT INFORMATION
MOBILE LEADED BARRIER WITH 36" X 24" WINDOW

