





Equipped with

- Leg/foot support joint with electric motor
- 2 x short side rails (seat section)
- 2 x side rails (back section)
- Integrated leg plate
- Removable rinsing basin
- Mobility (castors Ø 2")
- Hand and foot control

With memory function for sixteen positions, pelvic elevation, backrest, height and leg/foot support joint are electrically adjustable. All adjustments and memory functions can be operated via hand or foot control. Safe working load: 660 lbs.


Select the color



atlantic blue
10108470



shale grey
10108480



titan grey
10108490



milky coffee
10108500



lavender purple
10108510



stone grey
10108520

Select the leg support system

Foot supports
with integrated handle (pair)
10108100

Leg supports, Göpel type,
with integrated handle (pair)
10108110



Select leg section, clip-on

atlantic blue
10108410

shale grey
10108420

titan grey
10108430

milky coffee
10108440

lavender purple
10108450

stone grey
10108460

