


# MODEL DH4-8GD

## 8 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	HEIGHT
INTERIOR DIMENSIONS	15.75 IN	14.5 IN	46.75 IN	40 CM	36 CM	118 CM
EXTERIOR DIMENSIONS	19 IN	20.75 IN	71.5 IN	48 CM	52 CM	181 CM
CRATED DIMENSIONS	32 IN	32 IN	80 IN	81 CM	81 CM	203 CM

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

### SIDE VIEW



### FRONT VIEW



# SO-LOW



**NET WEIGHT : 230 LBS.**  
**SHIPPING WEIGHT: 300 LBS.**

### SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY ( REQUIRES FORKLIFT OR PALLET JACK )  
DOCK-TO-DOCK DELIVERY WITH LIFTGATE ( REQUIRES FORKLIFT OR PALLET JACK )  
INSIDE DELIVERY ( WHITE GLOVE DELIVERY )