

## MODEL DH4-8GD

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

<b>Capacity [ Cu. Ft. ]</b>	<b>6.17 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>175 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>
<b>Temperature Range [ °F ]</b>	<b>36°F to 46°F</b>
Dimensions	Width x Depth x Height
Exterior [ Inches ]	19" x 20.75" x 71.5"
Interior [ Inches ]	15.75" x 14.5" x 46.75"
Exterior [ Centimeters ]	48 x 52 x 181
Interior [ Centimeters ]	40 x 36 x 118
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R600a
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter ( Top Back Wall )
Casters	4
Storage	5 Adjustable Shelves
Shipping Weight	300 lbs.

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   2.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

#### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 3 Years Parts & Labor.  
Additional 4 Years on Compressor Parts, No Labor.