

## MODEL DH4-74SD

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

<b>Capacity [ Cu. Ft. ]</b>	<b>72 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>2095 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>
<b>Temperature Range [ °F ]</b>	<b>36°F to 46°F</b>
Dimensions	Width x Depth x Height
Exterior [ Inches ]	75" x 32" x 81"
Interior [ Inches ]	72" x 27" x 55"
Exterior [ Centimeters ]	190 x 81 x 206
Interior [ Centimeters ]	183 x 69 x 140
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Doors	Three Doors, Keyed Lock
Door Material	Insulated Steel
Interior Light	Adjustable LED
Porthole	1/2" Diameter ( Top Back Wall )
Casters	4, Locking
Storage	15 Adjustable Shelves
Shipping Weight	720 lbs.

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   7.2A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

#### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 7 Years Parts, Labor & Compressor Part Warranty.