

MODEL DH4-23SD

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	23 Cu. Ft.
Capacity [Liters]	650 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	36°F to 46°F

Dimensions	Width x Depth x Height
Exterior [Inches]	27" x 32" x 81"
Interior [Inches]	24" x 28.5" x 55"
Exterior [Centimeters]	68 x 81 x 206
Interior [Centimeters]	61 x 72 x 140

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Interior Light	Adjustable LED
Porthole	1/2" Diameter (Top Back Wall)
Casters	4, Locking
Storage	5 Adjustable Shelves
Shipping Weight	440 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 3.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 7 Years Parts, Labor & Compressor Part Warranty.