


MODEL DH4-27GD-SS

27 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


| | WIDTH | DEPTH | HEIGHT | WIDTH | DEPTH | HEIGHT |
|---------------------|-------|----------|--------|--------|--------|--------|
| INTERIOR DIMENSIONS | 27 IN | 28.5 IN | 55 IN | 69 CM | 72 CM | 140 CM |
| EXTERIOR DIMENSIONS | 30 IN | 33.75 IN | 81 IN | 76 CM | 85 CM | 206 CM |
| CRATED DIMENSIONS | 46 IN | 43 IN | 86 IN | 116 CM | 109 CM | 218 CM |

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

SIDE VIEW



FRONT VIEW



SO-LOW



NET WEIGHT : 385 LBS.
SHIPPING WEIGHT: 440 LBS.

SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY (REQUIRES FORKLIFT OR PALLET JACK)
DOCK-TO-DOCK DELIVERY WITH LIFTGATE (REQUIRES FORKLIFT OR PALLET JACK)
INSIDE DELIVERY (WHITE GLOVE DELIVERY)