

## MODEL DH4-23SD-SS

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

|                                 |                     |
|---------------------------------|---------------------|
| <b>Capacity [ Cu. Ft. ]</b>     | <b>23 Cu. Ft.</b>   |
| <b>Capacity [ Liters ]</b>      | <b>651 Liters</b>   |
| <b>Temperature Range [ °C ]</b> | <b>2°C to 8°C</b>   |
| <b>Temperature Range [ °F ]</b> | <b>36°F to 46°F</b> |

|                          |                        |
|--------------------------|------------------------|
| Dimensions               | Width x Depth x Height |
| Exterior [ Inches ]      | 27.25" x 33.75" x 81"  |
| Interior [ Inches ]      | 24" x 28.5" x 55"      |
| Exterior [ Centimeters ] | 69 x 85 x 206          |
| Interior [ Centimeters ] | 61 x 72 x 140          |

|                 |                           |
|-----------------|---------------------------|
| Compressor      | Hermetic Compressor       |
| Air Circulation | Forced Air                |
| Defrost Method  | Cycle                     |
| Refrigerant     | R290                      |
| Exterior        | Stainless Steel           |
| Interior        | Stainless Steel           |
| Insulation      | Polyurethane Foam         |
| Lockable Door   | 1 Door, Keyed Lock        |
| Door Material   | Insulated Stainless Steel |
| Interior Light  | Adjustable LED            |
| Porthole        | 1/2" Diameter             |
| Casters         | 4, Locking                |
| Storage         | 3 Adjustable Shelves      |
| Shipping Weight | 440 lbs.                  |

#### Electrical

|                         |                                |
|-------------------------|--------------------------------|
| Voltage:                | 115 Volts 60 Hertz 1 Phase     |
| Amperage Line   Running | 15 Amp Dedicated   3.0A        |
| NEMA Configuration      | NEMA 5-15                      |
| Certification(s)        | UL Listed<br>Energy Star Rated |

#### Controller

|                       |                         |
|-----------------------|-------------------------|
| Display Type          | Microprocessor, Digital |
| Alarms                | High/Low Temp           |
| Min / Max Temperature | Display and Reset       |
| Alarm Relay           | Dry Contacts            |

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 7 Years Parts, Labor & Compressor Part Warranty.