

# SO-LOW

ENVIRONMENTAL EQUIPMENT COMPANY

SETTING TEMPERATURE STANDARDS IN EVERY LAB

EST. 1959



## ULTRA LOW FREEZER & LAB REFRIGERATION EQUIPMENT PRODUCT CATALOG

MEDICAL RESEARCH | INDUSTRIAL USE | COLD STORAGE | LAB TESTING | LIFE SCIENCES



MEDICAL RESEARCH | LAB TESTING | VACCINE STORAGE | TEMPERATURE CRITICAL STORAGE

### Laboratory Quality Design and Operation

So-Low Glass Door Refrigerators maintain temperatures between 2°C to 8°C, and are primarily used for vaccines or storage of temperature critical items. Units come standard with High/Low alarms. Power failure and alarm dry contact relay are included.

### Temperature Critical Vaccine Storage

So-low lab refrigerators meet CDC recommendations for vaccine storage.

### UL Listed

UL listed to Standard 61010-1. Energy Star Rated.



MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	DEFROST TYPE	UNIT WEIGHT
--------------	------------	--------	-------------------	------------------------------------	------------------------------------	--------------	-------------

### GLASS DOOR UPRIGHT

DH4-8GD	6.17	175	+2°C to +8°C	15.75" x 15.5" x 46.75"	19" x 20.75" x 71.5	CYCLE	250 LBS
DH4-12GD	9.19	260	+2°C to +8°C	22" x 16.5" x 43.75"	25" x 22.75" x 63.25"	CYCLE	350 LBS
DH4-15GD	15.4	436	+2°C to +8°C	25.5" x 19.25" x 54.25"	29.5" x 26" x 78.5"	CYCLE	400 LBS
DH4-23GD	23	650	+2°C to +8°C	24" x 28.5" x 55"	27" x 32" x 81"	CYCLE	400 LBS
DH4-27GD	27	764	+2°C to +8°C	27" x 28.5" x 55"	30" x 32" x 81"	CYCLE	400 LBS
DH4-38GD	38	1060	+2°C to +8°C	43" x 22.25" x 60.5"	47.25" x 29.25" x 78.5"	CYCLE	450 LBS
DH4-49GD	49	1387	+2°C to +8°C	49" x 28.5" x 55"	52" x 32" x 81"	CYCLE	560 LBS
DH4-74GD	72	2095	+2°C to +8°C	72" x 27" x 61"	75" x 32" x 81"	CYCLE	700 LBS

### SLIDING GLASS DOOR UPRIGHT

DH4-30SGD	38	1060	+2°C to +8°C	43" x 23.25" x 60.5"	47.25" x 29.25" x 78.5"	CYCLE	450 LBS
DH4-38SGD	38	1076	+2°C to +8°C	40.5" x 27" x 55"	43.5" x 32" x 81"	CYCLE	500 LBS
DH4-45SGD	45	1286	+2°C to +8°C	49" x 27" x 55"	52" x 32" x 81"	CYCLE	560 LBS
DH4-66SGD	66	1868	+2°C to +8°C	72" x 27" x 55"	75" x 32" x 81"	CYCLE	800 LBS

### STAINLESS STEEL GLASS DOOR UPRIGHT

DH4-23GD-SS	23	651	+2°C to +8°C	24" x 28.5" x 55"	27.25" x 33.75" x 81"	CYCLE	390 LBS
DH4-27GD-SS	27	764	+2°C to +8°C	27" x 28.5" x 55"	30" x 33.75" x 81"	CYCLE	400 LBS
DH4-49GD-SS	49	1387	+2°C to +8°C	49" x 28.5" x 55"	52" x 33.75" x 81"	CYCLE	560 LBS
DH4-72GD-SS	72	2095	+2°C to +8°C	72" x 27" x 61"	75" x 33.75" x 81"	CYCLE	800 LBS