

SO-LOW

ENVIRONMENTAL EQUIPMENT COMPANY

SETTING TEMPERATURE STANDARDS IN EVERY LAB

EST. 1959



ULTRA LOW FREEZER & LAB REFRIGERATION EQUIPMENT PRODUCT CATALOG

MEDICAL RESEARCH | INDUSTRIAL USE | COLD STORAGE | LAB TESTING | LIFE SCIENCES



MEDICAL RESEARCH | LAB TESTING | VACCINE STORAGE | TEMPERATURE CRITICAL STORAGE

Laboratory Quality Design and Operation

So-Low Solid Door Refrigerators maintain temperatures between 2°C to 8°C, and are primarily used for vaccines or storage of temperature critical items. Units come standard with High/Low alarms. Power failure and alarm dry contact relay are included.

Photosensitive Protection

Solid doors block ambient light, allowing for safe storage of light sensitive products. Interior LED light can be enabled or disabled.

Temperature Critical Vaccine Storage

So-low lab refrigerators meet CDC recommendations for vaccine storage.

Additional Storage Safety

Temperature recorders, phone alarms, other safety options are available.

UL Listed

UL listed to Standard 61010-1. Energy Star Rated.



MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	DEFROST TYPE	UNIT WEIGHT
--------------	------------	--------	-------------------	---------------------------------	---------------------------------	--------------	-------------

SOLID DOOR UPRIGHT

DH4-23SD	23	650	+2°C to +8°C	24" x 28.5" x 55"	27" x 34" x 81"	CYCLE	390 LBS
DH4-27SD	27	764	+2°C to +8°C	27" x 28.5" x 55"	30" x 33.75" x 81"	CYCLE	400 LBS
DH4-49SD	49	1387	+2°C to +8°C	49" x 28.5" x 55"	52" x 33.75" x 81"	CYCLE	560 LBS
DH4-74SD	74	2095	+2°C to +8°C	72" x 27" x 61"	75" x 33.75" x 81"	CYCLE	800 LBS

STAINLESS STEEL SOLID DOOR UPRIGHT

DH4-23SD-SS	23	650	+2°C to +8°C	24" x 28.5" x 55"	27.25" x 28.5" x 55"	CYCLE	390 LBS
DH4-27SD-SS	27	764	+2°C to +8°C	27" x 28.5" x 55"	30" x 33.75" x 81"	CYCLE	400 LBS
DH4-49SD-SS	49	1387	+2°C to +8°C	49" x 28.5" x 55"	52" x 33.75" x 81"	CYCLE	560 LBS
DH4-72SD-SS	72	2095	+2°C to +8°C	72" x 27" x 61"	75" x 33.75" x 81"	CYCLE	800 LBS