

# SO-LOW

ENVIRONMENTAL EQUIPMENT COMPANY

SETTING TEMPERATURE STANDARDS IN EVERY LAB

EST. 1959



## ULTRA LOW FREEZER & LAB REFRIGERATION EQUIPMENT PRODUCT CATALOG

MEDICAL RESEARCH | INDUSTRIAL USE | COLD STORAGE | LAB TESTING | LIFE SCIENCES



**MEDICAL RESEARCH | LAB TESTING | VACCINE STORAGE | TEMPERATURE CRITICAL STORAGE**

### Touchscreen Control Display

The touchscreen series has an advanced PID control system, and a 4.3" color touch screen operating interface with embedded SCADA (supervisory control and data acquisition) and LAN features. These features include data logging, alarm management, profile capability, security options, operator audit trail, USB file transfer; as well as email/text alerts for alarm state.

### Temperature Critical Vaccine Storage


Units maintain temperatures between 2°C to 8°C, and meet CDC recommendations for vaccine storage.

### Local Network Remote Access

Features VNC remote access, for LAN remote control monitoring.

### UL Listed

UL listed to Standard 61010-1.



MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	DEFROST TYPE	UNIT WEIGHT
--------------	------------	--------	-------------------	------------------------------------	------------------------------------	--------------	-------------

#### GLASS DOOR UPRIGHT

DHP4-23SD	23	650	+2°C to +8°C	24" x 28.5" x 55"	27" x 34" x 81"	CYCLE	390 LBS
DHP4-27SD	27	764	+2°C to +8°C	27" x 28.5" x 55"	30" x 33.75" x 81"	CYCLE	400 LBS
DHP4-49GD	49	1387	+2°C to +8°C	49" x 28.5" x 55"	52" x 33.75" x 81"	CYCLE	560 LBS
DHP4-74GD	74	2095	+2°C to +8°C	72" x 27" x 61"	75" x 33.75" x 81"	CYCLE	700 LBS

#### SLIDING GLASS DOOR UPRIGHT

DHP4-38SGD	38	1076	+2°C to +8°C	40.5" x 27" x 55"	43.5" x 32" x 81"	CYCLE	500 LBS
DHP4-45SGD	45	1286	+2°C to +8°C	49" x 27" x 55"	52" x 32" x 81"	CYCLE	560 LBS
DHP4-66SGD	66	1868	+2°C to +8°C	72" x 27" x 55"	75" x 32" x 81"	CYCLE	800 LBS

#### SOLID DOOR UPRIGHT

DHP4-23GD	23	650	+2°C to +8°C	24" x 28.5" x 55"	27" x 34" x 81"	CYCLE	390 LBS
DHP4-27GD	27	764	+2°C to +8°C	27" x 28.5" x 55"	30" x 33.75" x 81"	CYCLE	400 LBS
DHP4-49SD	49	1387	+2°C to +8°C	47.5" x 28.5" x 55"	52" x 33.75" x 81"	CYCLE	560 LBS
DHP4-74SD	74	2095	+2°C to +8°C	72" x 27" x 55"	75" x 33.75" x 81"	CYCLE	800 LBS