

# AWC12PB



68.13" x 23.5" x 23.88" (H x W x D)

Performance Series hospital blanket warmer with a self-closing glass door, stainless steel interior, digital thermostat, casters, and lock

### Highlights:

Full length glass door keeps complete interior in view for easy inventory check

Factory installed lock allows you to restrict access to staff

Antimicrobial handle to help reduce germ spread

### Product Features:

Performance Series Blanket Warmer	Designed to keep blankets warm for patient care in hospitals, doctors' offices, physical therapy centers, and similar settings
Pre-set to 130°F	Standard temperature for keeping blankets warm for patient care (NOTE: parameters are locked at a 130°F maximum, but can be adjusted up to 160°F. Contact the Accucold technical support team for instructions)
Microprocessor control panel	Digital controls with a large LED readout in Celsius or Fahrenheit for intuitive temperature management
Antimicrobial silver ion-coated handle	Powder-coated handle with naturally occurring silver ions to help reduce the spread of germs
Double pane tempered glass door	Heat-safe glass door offers a full view of stored contents
Self-closing door	Door will automatically close when left ajar, preventing items from accidentally being left exposed to ambient temperatures
Stainless steel interior	Interior walls are constructed from polished 304 grade stainless steel to ensure lasting durability and easy cleaning
Adjustable stainless steel shelving	Heavy duty capacity (up to 44 lbs.) with a perforated design to help ensure proper heat distribution
Hospital grade casters	7" high rubberized casters (2 locking) ensure easy mobility as needed

## AWC12PB Specifications:

Overview	
Height of Cabinet	68.13" (173 cm)
Width	23.5" (60 cm)
Depth	23.88" (61 cm)
Depth with Handle	25.38" (64 cm)
Depth with door at 90°	45.63" (116 cm)
Capacity	12.54 cu.ft. (355 L)
Door	Glass
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	3.4
Voltage/Frequency	115 V AC/60 Hz
Weight	140.0 lbs. (64 kg)
Parts & Labor Warranty	2 Years
UPC	761101150994
Warming Cabinet	
Door Swing	RHD
Reversible	Factory Reversible
Adjustable Shelves	Yes
Shelf Type	Steel
Shelf Qty	5
Thermostat Type	Digital
Interior Light	Yes
Interior Height	56.0" (142 cm)
Interior Width	19.75" (50 cm)
Interior Depth	18.0" (46 cm)
Wheels	4
Temperature Range	90 to 130°F; To adjust up to 160°F, contact our technical support team



# AWC12PBLHD



68.13" x 23.5" x 23.88" (H x W x D)

Performance Series hospital blanket warmer with a self-closing glass door, stainless steel interior, digital thermostat, casters, and lock

### Highlights:

Full length glass door keeps complete interior in view for easy inventory check

Factory installed lock allows you to restrict access to staff

Antimicrobial handle to help reduce germ spread

### Product Features:

Performance Series Blanket Warmer	Designed to keep blankets warm for patient care in hospitals, doctors' offices, physical therapy centers, and similar settings
Pre-set to 130°F	Standard temperature for keeping blankets warm for patient care (NOTE: parameters are locked at a 130°F maximum, but can be adjusted up to 160°F. Contact the Accucold technical support team for instructions)
Microprocessor control panel	Digital controls with a large LED readout in Celsius or Fahrenheit for intuitive temperature management
Antimicrobial silver ion-coated handle	Powder-coated handle with naturally occurring silver ions to help reduce the spread of germs
Double pane tempered glass door	Heat-safe glass door offers a full view of stored contents
Self-closing door	Door will automatically close when left ajar, preventing items from accidentally being left exposed to ambient temperatures
Stainless steel interior	Interior walls are constructed from polished 304 grade stainless steel to ensure lasting durability and easy cleaning
Adjustable stainless steel shelving	Heavy duty capacity (up to 44 lbs.) with a perforated design to help ensure proper heat distribution
Hospital grade casters	7" high rubberized casters (2 locking) ensure easy mobility as needed

## AWC12PBLHD Specifications:

Overview	
Height of Cabinet	68.13" (173 cm)
Width	23.5" (60 cm)
Depth	23.88" (61 cm)
Depth with Handle	25.38" (64 cm)
Depth with door at 90°	45.63" (116 cm)
Capacity	12.54 cu.ft. (355 L)
Door	Glass
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	3.4
Voltage/Frequency	115 V AC/60 Hz
Weight	140.0 lbs. (64 kg)
Parts & Labor Warranty	2 Years
UPC	761101155111
Warming Cabinet	
Door Swing	LHD
Reversible	No
Adjustable Shelves	Yes
Shelf Type	Steel
Shelf Qty	5
Thermostat Type	Digital
Interior Light	Yes
Interior Height	56.0" (142 cm)
Interior Width	19.75" (50 cm)
Interior Depth	18.0" (46 cm)
Wheels	4
Temperature Range	90 to 130°F; To adjust up to 160°F, contact our technical support team

