

QUICKSHIP LEAD APRON TRAVEL BAG

PRODUCT INFORMATION



PRODUCT INFORMATION

QUICKSHIP LEAD APRON TRAVEL BAG | QS-RA-BAG-BK



The Phillips Safety QuickShip Lead Apron Bag is the essential travel solution for radiologists, doctors, and technicians. Designed to preserve the integrity of your apron's internal lining, this bag eliminates the risk of cracks and creases caused by improper folding.

Key Features:

Crease Prevention: Includes a heavy-duty tube to roll your apron around, ensuring long-term shielding safety.

Easy Transport: Features both a reinforced grab handle and an adjustable shoulder strap for hands-free mobility.

QuickShip: In-stock and ready for immediate dispatch to meet your tight schedule.

Durable Build: Constructed from rugged, easy-to-clean Black Nylon.

Medical Grade: 100% Latex-Free.

Simple Storage: Roll your apron around the provided tube. Slide the rolled apron into the bag. Zip and go. Stop folding your protection. Maintain your apron's lifespan with a professional storage solution that's ready to ship when you are.

SPECIFICATION

PRODUCT TYPE: Travel Bag

MATERIAL: Nylon

COLOR: Black

FEATURES: Reinforced grab handle; Adjustable shoulder

strap; Heavy-duty tube; Easy-to-Clean; Latex-Free