

LEAD APRON TRAVEL BAG

PRODUCT INFORMATION

PRODUCT INFORMATION
LEAD APRON TRAVEL BAG | RA-BAG-CUST



The Phillips Safety Custom Lead Apron Bag is the essential travel solution for radiologists, doctors, and technicians. Designed to preserve the integrity of your apron’s internal lining, this bag eliminates the risk of cracks and creases caused by improper folding.

Key Features:

- Crease Prevention:** Includes a heavy-duty tube to roll your apron around, ensuring long-term shielding safety.
- Easy Transport:** Features both a reinforced grab handle and an adjustable shoulder strap for hands-free mobility.
- Durable Build:** Constructed from rugged, easy-to-clean Black Nylon.
- Medical Grade:** 100% Latex-Free.
- Customized Options:** You can choose the color/pattern, trim, and embroidery of your bag.

Simple Storage: Roll your apron around the provided tube. Slide the rolled apron into the bag. Zip and go. Stop folding your protection. Maintain your apron’s lifespan with a professional storage solution that’s ready to ship when you are.

SPECIFICATION

PRODUCT TYPE: Travel Bag	FEATURES: Reinforced grab handle; Adjustable shoulder strap; Heavy-duty tube; Easy-to-Clean; Latex-Free
MATERIAL: Nylon	CUSTOMIZED OPTIONS: Color/Pattern, Trim, and Embroidery