

AFS14PV



74.38" x 23.5" x 25.25" (H x W x D)

Performance series Pharma-Vac 9.53 cu. ft. upright frost-free all-freezer for vaccine storage

Highlights:

Designed and purpose built to meet demanding needs of life science, vaccine, laboratory and medical cold storage applications

Intelligent microprocessor digital temperature controller

True frost-free operation decreases user maintenance by preventing icy buildup

Product Features:

24" Wide Upright All-Freezer	24" Wide upright Performance Series vaccine all-freezer designed for healthcare settings
Performance Series freezer	Purpose-built all-freezer for vaccine storage in medical, pharmaceutical, and life science applications
Purpose-built design	Designed and purpose-built for pharmacy, medication, and vaccine applications to support meeting CDC/VFC/VFA guidelines
Spacious capacity	Over 9.5 cu.ft. of storage inside
Frost-free operation	No-frost operation for minimal user maintenance
+/-1°C Variation	The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period
Microprocessor control panel	Digital controls with a large LED readout in Celsius or Fahrenheit for intuitive temperature management
Adjustable temperature	Intuitive controls let you adjust the interior temperature between -30°C to -10°C
High/low temperature alarm	Visual and audible alarms alert if the interior temperature rises or falls out of set range

AFS14PV Specifications:

Overview	
Nominal Height	75.0" (191 cm)
Height of Cabinet	74.38" (189 cm)
Nominal Width	24.0" (61 cm)
Width	23.5" (60 cm)
Width with Door Open	23.75" (60 cm)
Nominal Depth	26.0" (66 cm)
Depth	25.25" (64 cm)
Depth with Handle	26.75" (68 cm)
Depth with door at 90°	49.38" (125 cm)
Capacity	9.53 cu.ft. (270 L)
Defrost Type	Frost-Free
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	4.3
Voltage/Frequency	115 V AC/60 Hz
Parts & Labor Warranty	2 Years
Compressor Warranty	5 Years
UPC	761101147185
Freezer	
Door Swing	RHD
Reversible	Yes
Interior Height	43.75" (111 cm)
Interior Width	15.25" (39 cm)
Interior Depth	16.5" (42 cm)
Shelf Type	Wire
Shelf Qty	6
Adjustable Shelves	Yes
Thermostat Type	Digital
Interior Light	Yes
Fan Type	Interior
Refrigerant Type	R290a
Refrigerant Amount	3.27 oz.
High Side PSI	280.0
Low Side PSI	130.0
Wheels	4

Compressor Step Height	10.0" (25 cm)
Compressor Step Width	15.25" (39 cm)
Compressor Step Depth	15.25" (39 cm)
Temperature Range	-30 to -10°C



AFS14PVLHD



74.38" x 23.5" x 25.25" (H x W x D)

Performance series Pharma-Vac 9.53 cu. ft. upright frost-free all-freezer for vaccine storage

Highlights:

Designed and purpose built to meet demanding needs of life science, vaccine, laboratory and medical cold storage applications

Intelligent microprocessor digital temperature controller

True frost-free operation decreases user maintenance by preventing icy buildup

Product Features:

24" Wide Upright All-Freezer	24" Wide upright Performance Series vaccine all-freezer designed for healthcare settings
Performance Series freezer	Purpose-built all-freezer for vaccine storage in medical, pharmaceutical, and life science applications
Purpose-built design	Designed and purpose-built for pharmacy, medication, and vaccine applications to support meeting CDC/VFC/VFA guidelines
Spacious capacity	Over 9.5 cu.ft. of storage inside
Frost-free operation	No-frost operation for minimal user maintenance
+/-1°C Variation	The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period
Microprocessor control panel	Digital controls with a large LED readout in Celsius or Fahrenheit for intuitive temperature management
Adjustable temperature	Intuitive controls let you adjust the interior temperature between -30°C to -10°C
High/low temperature alarm	Visual and audible alarms alert if the interior temperature rises or falls out of set range

AFS14PVLHD Specifications:

Overview	
Nominal Height	75.0" (191 cm)
Height of Cabinet	74.38" (189 cm)
Nominal Width	24.0" (61 cm)
Width	23.5" (60 cm)
Width with Door Open	23.75" (60 cm)
Nominal Depth	26.0" (66 cm)
Depth	25.25" (64 cm)
Depth with Handle	26.75" (68 cm)
Depth with door at 90°	49.38" (125 cm)
Capacity	9.53 cu.ft. (270 L)
Defrost Type	Frost-Free
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	4.3
Voltage/Frequency	115 V AC/60 Hz
Parts & Labor Warranty	2 Years
Compressor Warranty	5 Years
UPC	761101152219
Freezer	
Door Swing	LHD
Interior Height	43.75" (111 cm)
Interior Width	15.25" (39 cm)
Interior Depth	16.5" (42 cm)
Shelf Type	Wire
Shelf Qty	6
Adjustable Shelves	Yes
Thermostat Type	Digital
Interior Light	Yes
Fan Type	Interior
Refrigerant Type	R290a
Refrigerant Amount	3.27 oz.
High Side PSI	280.0
Low Side PSI	130.0
Wheels	4
Compressor Step Height	10.0" (25 cm)

Compressor Step Width	15.25" (39 cm)
Compressor Step Depth	15.25" (39 cm)
Temperature Range	-30 to -10°C

