

CAMERA COMPARISON CHART



MODEL	MegaVID WiFi Camera	BioVID 4K Camera & Monitor	BioVID 4K	BioVID HD Camera & Monitor	BioVID HD 1080+	MiniVID USB 3.0 / 6.3MP
MODEL NUMBER	MVC-12MP-WiFi	BVC-4K16-TVK3	BVC-4K16-CMT3	BVC-1080-TVK3	BVC-1080-CMT3	MVC-U5MP-EMTN
RESOLUTION	12MP	4K - 8MP	4K - 8MP	1080 - 2MP	1080 - 2MP	6.3MP
VIDEO OUTPUT	HDMI, USB, WIFI	HDMI, USB	HDMI, USB	HDMI, USB2	HDMI, USB2	USB3
OPERATING SYSTEM	PC, Mac, iOS, and Android	PC or MAC	PC or MAC	PC or MAC	PC or MAC	PC or MAC
MOUNTING OPTIONS	C-Mount	C-Mount	C-Mount	C-Mount	C-Mount	C-Mount, Eyetube
KEY FEATURES	High Definition 12MP resolution, Connects to TV, PC, MAC, Phone, and Tablet	HD imaging to 13.3" monitor. Save to flash drive.	Ultra High Definition 4K resolution. Connects to TV, and PC/MAC.	High Definition 1080p imaging to TV. Capture images & video to computer or SD card.	High Definition 1080p imaging to TV. Capture images & video to computer or SD card.	Sony sensor; 3.0 USB for super-fast "LIVE" streaming, automatic exposure, ToupView micro-image analysis software.
APPLICATIONS	High resolution live presentations, digital capture, teaching & client education	Highest resolution camera for universities, laboratories, and veterinary applications	Highest resolution camera for universities, laboratories, and veterinary applications	High resolution live presentations, digital capture, teaching, & client education	High resolution live presentations, digital capture, teaching, & client education	Digital capture for education, documentation, email