

# MEGAVID

## WI-FI CAMERA

HIGHEST  
RESOLUTION  
FASTEST  
FRAME RATE!

MORE WAYS TO CONNECT  
TO IMAGES AND PEOPLE.




HDMI



USB





WI-FI



ETHERNET

STREAM LIVE IMAGES TO ALL YOUR ANDROID AND IOS DEVICES  
AT THE SAME TIME!



## IOS AND ANDROID PHONES AND TABLETS · PC COMPUTERS · 4K TV'S AND PROJECTORS

The MegaVID Wi-Fi 12MP camera is an incredible addition to our camera line. Operating at 5GHz frequency for maximum speed (4000 x 3000 @ 30fps) the MegaVID can be mounted to nearly any microscope that supports trinocular C-mount and connects to any iOS or Android device through the LW Scientific app. Utilize the MegaVID Wi-Fi to its fullest capabilities using the included software disc to connect to any desktop, any PC via USB, and any 4K TV or Projector via HDMI. The MegaVID Wi-Fi has many applications but is best suited for teaching, training, and educating clients and students.

# SPECIFICATIONS

- 12 megapixel camera 4000 x 3000 @ 30fps
- 1/2.3" Sony IMX412 CMOS sensor with 1.55µm pixels
- 5GHz high-speed Wi-Fi with H.264/MPEG-4 video compression
- Real time auto and manual white balance adjustment
- 5vDC USB power supply

# FEATURES

- Live stream video onto any hand-held device, PC, or HDTV
- High-res imaging in BMP, JPG, and PNG formats
- LW Scientific App available on iOS and Android
- LW Scientific Software included for PC

# REQUIREMENTS

- iOS11 or higher
- Android 3.0 or higher
- Microsoft Windows 7 or higher (32 & 64 bit)
- Compatible C-Mount Trinocular pipe



# PACKAGE INCLUDES

- [1] MegaVID Wi-Fi camera
- [1] HDMI Cable
- [1] Ethernet Cable
- [1] Software Disc
- [1] Power Adapter



SKU MVC-12MP-WIFI

