

# MiniVID USB Camera

The MiniVID USB is our best-selling camera for microscopy. It can be mounted to nearly any brand or type of microscope, either into the eye-tube of a binocular microscope, or thread-mounted onto any trinocular C-mount. Choose from two digital models: USB3.0 or USB2.0. The MiniVID USB is the perfect addition to any microscope for training techs, educating clients, or documenting for patient files and medical device records.

## **MINIVID 6.3mp USB3.0 - Highest Quality for Professional Use**

Very fast frame rate for fast-action digital video. Works on any MAC or PC computer with USB3.0 port.



## MINIVID 6.3MP USB 3.0

### Key Features.

Sony sensor; 3.0 USB for super-fast "LIVE" streaming

### Image Sensor

1/1.8" Color Sony CMOS chip - 6.3mp

### Resolution

3072 X 2048, 2.4um pixel size

### Interface

USB 3.0

### System Requirement

Windows 98, 7, 8, XP, 10, Vista, Mac OS

### Software

ToupView micro-image analysis software  
Save digital images and videos

### Mounting Requirements

C-Mount, Eyetube (optical adapter included)

### Frame Size and Rate

30fps @ 3072 x 2048, 38fps @ 1536 x 1024

## MiniVID USB Camera

Model #	Description
MVC-U6MP-USB3	6.3mp USB3.0 Super-Fast Camera w/ Sony Sensor and software

