

Fluorescent and LED Ring Lights

Fluorescent and LED ring lights offer shadow-free cold illumination with super-long-life bulbs at a reduced cost. Long thumb screws allow for easy attachment to almost any brand of microscope. The newly-redesigned fluorescent ring light is the most affordable option, and the LED models produce the most intense, focused illumination with 10,000-hour white LED bulbs. The economical 48-bulb LED ring light features variable intensity control, and the new 60-bulb LED ring light is variable and programmable for segmented quadrants to cast shadows when necessary.



Fluorescent Ring Light

FLUORESCENT

10-watt, 6400K bulb
 110v (220v) AC input
 Spare bulb included
 Fits 33mm-69mm mount rings
 Boxed: 0.8 lb
 11" x 8.3" x 2.25"



48-Bulb Segmented Ring Light

SEGMENTED 48-BULB LED

Segmented quadrants for shadowing
 48 bright white focused LED bulbs
 110v / 220v auto switching
 1 amp fuse
 Adjustable dimmer
 Fits 30mm-60mm mount rings
 Boxed: 0.8 lb
 9" x 8" x 2"

Fluorescent and LED Ring Lights

Model #	Description
ILL-FLR7-1071	Fluorescent Ring Light 10w (non-variable)
ILL-LEDV-R483	Segmented Variable LED Ring Light - 48 bulb
ILL-LEDS-R603	Segmented Variable LED Ring Light - 60 bulb

SEGMENTED 60-BULB LED

Segmented quadrants for shadowing
 60 bright white focused LED bulbs
 110v / 220v auto switching
 1 amp fuse
 Adjustable dimmer
 Fits 30mm-60mm mount rings
 Boxed: 0.8 lb
 9" x 8" x 2"

60-Bulb Segmented Ring Light

