

INJURY SPECIFIC FIRST AID

POWERED BY

Modulator 



Emergency Response Kit


Emergency Response | Softpack

The Modulator Ansi Class B compliant Modulator Emergency Response Kit is designed for ultimate preparedness, offering well-organized, injury-specific first aid modules for burns, eye injuries, sprains & fractures, bleed control, general first aid + an AED/CPR prep kit. Housed in a durable, weather-resistant soft case, this kit is built to withstand the toughest conditions, making it essential for any environment. Portable and reliable, the Modulator Emergency Response Kit is your trusted partner in emergency situations—ready to save lives, anytime, anywhere.



Code: M610B
Modulator Emergency Response Bag Standard

Code: M610B - 1
Modulator Emergency Response Bag Pro
This kit includes an upgraded pro bleed control kit.



Kit Includes a Metal Hook with Adhesive Backing

Aero's durable metal hook secures the Emergency Response Kit beneath a first aid cabinet or mounts directly to the wall keeping critical first aid supplies organized and easily accessible in one centralized location.

Emergency Response | Softpack Refills



BURN MODULE REFILL

Code: M6001 Case Included

CODE	CONTENTS	QTY
ABD10	AeroBurn Burn Dressing 4in x 4in	2
AF100	AeroForm Conforming Bandage 4in x 4.3yds	1
AGNPF02	AeroGlove Nitrile Gloves Large	2
APD57550BS	AeroPad Low Adherent Dressing 2in x 3in	2
ASC13PUS	AeroSupplies Scissors Blue Lister 5.2in	1
AB0905B	AeroBurn Burn Gel Packet 3.5g	25
ABS6090	AeroBurn PE Burn Sheet 23 1/2 x 35 1/2in	1
ABS2020	AeroBurn PE Burn Sheet 8in x 8in	1
ABS1010	AeroBurn PE Burn Sheet 4in x 4in	1
6001BG	Zippered Bag	1



EYE MODULE REFILL

Code: M6002 Case Included

CODE	CONTENTS	QTY
AEP15	AeroPad Eye Pads 2in x 3in	4
AF25	AeroForm Conforming Bandage 1in x 4.4yds	1
AF75	AeroForm Conforming Bandage 3in x 4.4yds	1
AGNPF02	AeroGlove Nitrile Gloves Large	2
AEM6PSUS	AeroSupplies Eye Magnet Plastic	1
AW1000US	AeroWash Eyewash .5oz	8
AP312	AeroTape Microporous Paper Tape .5in x 5.5yd	1
6002BG	Eye Wound Zippered Bag	1



GENERAL FIRST AID MODULE

Code: M6006 Case Included

CODE	CONTENTS	QTY
AFP60125	AeroPlast Premium Fabric X-Wide Strip 1in x 3in	25
AP10125	AeroPlast Plastic Strip 0.7in x 2.8in	25
AFP312	AeroPlast Premium Fabric Fingertip & Knuckle Dressings	12
AF50	AeroForm Conforming Bandage 2in x 4.3yds	2
AF9	AeroSupplies Forceps Steel 3.5in	1
ASP100	AeroProbe Splinter Probes 1.4in	5
ASC13PUS	AeroSupplies Scissors Blue Lister 5.2in	1
AW8001	AeroWipe Alcohol Swab	30
AW7001	Alcohol Wipe 1.25in x 1.25in	10
AGNPF02	AeroGlove Nitrile Gloves Large 1 Pair	1
AT25	AeroPlast Transparent Tape 1in x 5.5yd	1
APD57550BS	AeroPad Low Adherent Dressing 2in x 3in	1
APG1050BSUS	AeroPad Gauze Pad 4in x 4in	2
APG7550BS	AeroPad Gauze Pad 3in x 3in	4
APG550BSUS	AeroPad Gauze Pad 2in x 2in	2
AGG001	AeroGuide First Aid Guide	1
ASTA05US	AeroSoothe Triple Antibiotic 0.5g	25
ACG10US	AeroCleanse Hand Sanitizer 0.9g	20
ACD1223S	AeroWound Combine Dressing 5in x 9in	1
AWF002	AeroWound Finger Dressing	1
M6006_BAG	General First Aid Zippered Bag	1

This sturdy metal wall bracket keeps your response kit visible and ready for emergencies.



Item: Metal Wall Bracket

Code: MWB01

Emergency Response | Softpack Refills



SPRAINS & FRACTURES MODULE

Code: M6007 Case Included

CODE	CONTENTS	QTY
ATB213	AeroRescue Thermal Blanket 55in x 84in	1
ABN10	AeroBand Triangular Bandage	2
ASP12	AeroSupplies Safety Pins Medium	12
AGNPF02	AeroGlove Nitrile Gloves Large 1 Pair	1
AFH100	AeroForm Conforming Heavyweight Bandage 4in x 4.3yds	1
AFH50	AeroForm Conforming Heavyweight Bandage 2in x 4.3yds	1
All3100US	AeroPlast Instant Ice Pack 4in x 5in	2
AW8001	AeroWipe Antiseptic Wipe	10
AR1060R	AeroSplint Aluminium Foam Splint 4.25in x 24in Rolled	1
M6007_BAG	Sprains & Fractures Zippered Bag	1



STANDARD BLEED CONTROL MODULE

Code: M6008 Case Included

CODE	CONTENTS	QTY
AF100	AeroForm Conforming Bandage 4in x 4.3yds	1
ACD1223S	AeroWound Combine Dressing 5in x 9in	3
AW8001	AeroWipe Antiseptic Wipe 5in x 7in	10
RST103	RapidStop Arterial Tourniquet	1
RSC104	Permanent Marker Black	1
AGNPF02	AeroGlove Nitrile Gloves Large 1 Pair	2
ATB213	AeroRescue Thermal Blanket 55in x 84in	1
AWD100S	AEROWOUND Pressure Dressing 4in x 4in	1
AWD75SS	AEROWOUND Pressure Dressing 3in x 3in	1
AWD50S	AEROWOUND Pressure Dressing 2in x 2in	1
M6008_BAG	Bleed Control Zippered Bag	1



PRO BLEED CONTROL MODULE

Code: M6009 Case Included

CODE	CONTENTS	QTY
AGNPF02	AeroGlove Nitrile Gloves Large 1 Pair	2
ATB213	AeroRescue Thermal Blanket 55in x 84in	1
RST103	RapidStop Arterial Tourniquet	1
RSC104	Permanent Marker Black	1
RSC100	Hemorrhage Control Bandage 4in	1
RSC103	Chest Seal Dressing Vent 2/pack	2
RSC102	Hemostatic Gauze 3in x 4yd	1
USP15	Trauma Shears 5.5in	1
ACD1223S	AeroWound Combine Dressing 5in x 9in	3
AF100	AeroForm Conforming Bandage 4in x 4.3yds	1
AWD100S	AEROWOUND Pressure Dressing 4in x 4in	1
AWD75SS	AEROWOUND Pressure Dressing 3in x 3in	1
AWD50S	AEROWOUND Pressure Dressing 2in x 2in	1
AW8001	AeroWipe Antiseptic Wipe	10
M6008_BAG	Bleed Control Zippered Bag	1



AED / CPR PREP MODULE

Code: M6010 Case Included

CODE	CONTENTS	QTY
AM01BUS	AeroShield CPR Mask & Gloves	1
ARZ001	AeroSupplies Double Sided Surgical Razor	1
AGNPF02	AeroGlove Nitrile Gloves Large 1 Pair	2
AHBB24US	AeroHazard Biohazard Bag 24in x 24in	1
AW8001	AeroWipe Antiseptic Wipe 5in x 7in	1
USP15	AeroSupplies Universal Shears Plastic Handle 5.9in	1
AHPT	AeroHazard Absorbant Towel 10in x 13in	1
M6010_BAG	Yellow Zippered Bag	1



Emergency Response | Rugged Hard Case

The Rugged Hard Case Kits are designed to equip rescuers with the critical tools needed to handle everything from minor injuries to severe trauma, including cardiac arrest. Kits are available with the RapidStop® Tourniquet alone or can be upgraded to include a comprehensive Bleed Control Module that includes the RapidStop® Tourniquet along with a full range of bleeding control supplies. Both options are available with multiple AED configuration choices as detailed in the tables below.



Rugged Hard Case Kits With RapidStop® Tourniquet Only

CODE	CONFIGURATION
M600-1	Modulator Trauma Kit with Heartsine 350P - XL Rugged Hard Case
M600-2	Modulator Trauma Kit with Heartsine 360P - XL Rugged Hard Case
M600-3	Modulator Trauma Kit with Heartsine 450P - XL Rugged Hard Case
M600	Modulator Trauma Kit without AED - XL Rugged Hard Case



Rugged Hard Case Kits With Bleed Control Module

CODE	CONFIGURATION
M600-4	Modulator Trauma Kit with Heartsine 350P & Bleed Control - XL Rugged Hard Case
M600-5	Modulator Trauma Kit with Heartsine 360P & Bleed Control - XL Rugged Hard Case
M600-6	Modulator Trauma Kit with Heartsine 450P & Bleed Control - XL Rugged Hard Case
M600-7	Modulator Trauma Kit with Bleed Control without AED - XL Rugged Hard Case

Trauma & First Aid Kits



RapidStop Modulator Trauma Kit ANSI Class B – Rugged Case

The M610 is an exclusive Modulator Trauma & First Aid Kit meticulously crafted to equip rescuers with the necessary tools for addressing both major and minor injuries. At the heart of this comprehensive kit is the proprietary RapidStop® Tourniquet, which employs advanced technology for dependable, rapid occlusion. The M610 includes a rugged, weather-resistant carrying case, ensuring your first aid tools and supplies remain protected in extreme weather conditions and remote locations.

Code: M610



Modulator Trauma Kit with Bleed Control – Rugged Hard Case

While the M610 featured the revolutionary RapidStop® Tourniquet, the M610-1 goes a step further by incorporating our comprehensive **Bleed Control module**. This module, along with the trusted RapidStop® Tourniquet, brings together advanced technology for effective and rapid management of severe bleeding, ensuring dependable occlusion in critical scenarios.

Code: M610-1

All Rugged Hard Cases are waterproof for added protection.

All Modulator Kits Meet & Exceed:
ANSI / ISEA Z308.1-2021 Class B Type III Container
OSHA 1910.151(b)
OSHA 1910.266 Logging First Aid Kit
Cal/OSHA - Title 8, Section 1512



BURN MODULE REFILL
Code: M6001 Case Included



EYE MODULE REFILL
Code: M6002 Case Included



AED / CPR PREP KIT
Code: CA-PK2US Case Included

CODE	CONTENTS	QTY
ABD10	AeroBurn Burn Dressing 4in x 4in	2
AF100	AeroForm Conforming Bandage 4in x 4.3yds	1
AGNPF02	AeroGlove Nitrile Gloves Large	2
APD57550BS	AeroPad Low Adherent Dressing 2in x 3in	2
ASC13PUS	AeroSupplies Scissors Blue Lister 5.2in	1
AB0905B	AeroBurn Burn Gel Packet 3.5g	25
ABS6090	AeroBurn PE Burn Sheet 23 1/2 x 35 1/2in	1
ABS2020	AeroBurn PE Burn Sheet 8in x 8in	1
ABS1010	AeroBurn PE Burn Sheet 4in x 4in	1
6001BG	Zippered Bag	1

CODE	CONTENTS	QTY
AEP1S	AeroPad Eye Pads 2in x 3in	4
AF25	AeroForm Conforming Bandage 1in x 4.4yds	1
AF75	AeroForm Conforming Bandage 3in x 4.4yds	1
AGNPF02	AeroGlove Nitrile Gloves Large	2
AEM6PSUS	AeroSupplies Eye Magnet Plastic	1
AW1000US	AeroWash Eyewash .5oz	8
AP312	AeroTape Microporous Paper Tape .5in x 5.5yd	1
6002BG	Eye Wound Zippered Bag	1

CODE	CONTENTS	QTY
AM01BUS	AeroShield CPR Mask & Gloves	1
ARZ001	AeroSupplies Double Sided Surgical Razor	1
AGNPF02	AeroGlove Nitrile Gloves Large	4
AHBB24US	AeroHazard Biohazard Bag 24in x 24in	1
AW8001	AeroWipe Antiseptic Wipe 5in x 7in	1
USP15	AeroSupplies Universal Shears Plastic Handle 5.9in	1
AHPT	AeroHazard Absorbant Towel 10in x 13in	1
6005BG	AED / CPR Zippered Bag with Lanyard Clip	1

Emergency Response | Rugged Hard Case Refills



MINOR WOUND MODULE REFILL
Code: M6003 Case Included

CODE	CONTENTS	QTY
AFP60125	AeroPlast Premium Fabric X-Wide Strip 3in x 1in	25
AP10125	AeroPlast Plastic Strip 2.8 x 0.7in	25
AFP312	AeroPlast Premium Fabric Fingertip & Knuckle Dressings	12
AF50	AeroForm Conforming Bandage 2in x 4.3yds	2
AF9	AeroSupplies Forceps Steel 3.5in	1
ASP100	AeroProbe Splinter Probes 1.4in	5
AW8001	AeroWipe Antiseptic Wipe 5in x 7in	30
AW7001	Alcohol Swabs 1.25in x 1.25in	10
AGNPF02	AeroGlove Nitrile Gloves Large	2
AT25	AeroPlast Transparent Tape 1inx5.5yd	1
APD57550BS	AeroPad Low Adherent Dressing 2in x 3in	1
APG1050BSUS	AeroPad Gauze Pad 4in x 4in	2
APG7550BS	AeroPad Gauze Pad 3in x 3in	4
APG550BSUS	AeroPad Gauze Pad 2in x 2in	2
AGG001	AeroGuide First Aid Guide	1
ASTA05US	AeroSoothe Triple Antibiotic 0.5g	25
ACG10US	AeroCleanse Hand Sanitizer 0.9g	20
ACD1223S	AeroWound Combine Dressing 5in x 9in	1
AWF002	AeroWound Finger Dressing	1
6003BG	Minor Wound Zippered Bag	1



SERIOUS WOUND MODULE REFILL
Code: M6004 Case Included

CODE	CONTENTS	QTY
ATB213	AeroRescue Thermal Blanket 55in x 84in	1
ABN10	AeroBand Triangular Band. Non-Woven	2
ASP12	AeroSupplies Safety Pins Medium	12
AGNPF02	AeroGlove Nitrile Gloves Large	2
AF100	AeroForm Conforming Bandage 4in x 4.3yds	1
AFH100	AeroForm Conforming Heavyweight Bandage 4in x 4.3yds	1
AFH50	AeroForm Conforming Heavyweight Bandage 2in x 4.3yds	1
All3100US	AeroPlast Instant Ice Pack 4in x 5in	2
AW8001	AeroWipe Antiseptic Wipe 5in x 7in	20
ACD1223S	AeroWound Combine Dressing 5in x 9in	2
AWD15S	AeroWound Bandage Wound Dressing Large	1
6004BG	Serious Wound Zippered Bag	1



BLEED CONTROL MODULE REFILL*
Code: M6005 Case Included

CODE	CONTENTS	QTY
AGNPF02	AeroGlove Large Nitrile Powder-Free Gloves Pair/2	2
ATB213	AeroRescue Thermal Blanket 55in x 84in	1
RSC104	Permanent Marker Black	1
RSC100	Hemorrhage Control Bandage 4in	1
RSC103	Chest Seal Dressing Vent 2/pack	1
RSC102	Hemostatic Gauze 3in x 4yd	1
USP15	Trauma Shears 5.5in	1
RST103	RapidStop Military Tourniquet	1

*The Bleed Control Module is available as an additional purchase and is not included in the standard Modulator Trauma & First Aid Kit contents.

Emergency Response | Individual Items



Item: RapidStop® Tourniquet
Type: Responder
Color: Red
Code: RST103



Item: Padded Splint
Code: AR1060R
Size: 24"L x 4.25"H



Item: Metal Wall Bracket
Code: MWB01

Modulator 
injury specific first aid

YOUR ENTIRE EMERGENCY SAFETY NET In One Ultra Rugged Backpack

The **Modulator Trauma Backpack** by Aero offers ultimate preparedness with well-organized, injury-specific first aid modules for burns, eye injuries, sprains & fractures, bleed control, and general first aid, plus an AED/CPR prep kit. Contained within a rugged, weather-resistant backpack, this kit is built to endure the most challenging conditions, making it an essential tool for any environment.



Emergency Response | Backpack



Internal elastic straps securely organize and hold all your additional medical equipment.

A clear pouch allows extra storage with easy visibility.

Includes an extra empty module case to add your own special supplies.

Emergency Response | Backpack



The clear window allows a quick visual check of your AED - ensuring it's in the backpack and in a rescue ready state.

Choose Your Configuration

Select the level of preparedness that meets your needs. Our Standard backpack provides a comprehensive, ANSI Class B compliant foundation. The Pro version upgrades to the Pro Bleed Kit, which includes hemostatic and chest seal dressings for advanced hemorrhage control. For ultimate response capability, the Pro & Airway backpack adds a complete airway and resuscitation module to the Pro configuration.

CODE	CONFIGURATION
M600B-1	Modulator Trauma Backpack, without AED, Standard
M600B-2	Modulator Trauma Backpack, without AED, Pro
M600B-3	Modulator Trauma Backpack, without AED, Pro & Airway

A robust side pouch securely holds things like an oxygen bottle or a small fire extinguisher.



M600B-3

Optional Airway & Resuscitation Module provides a complete airway and resuscitation tool kit.



Features a dedicated pouch that fits most AEDs, including HeartSine, Philips, LIFEPAK CR2, ZOLL, and Defibtech



Adjustable straps and waste band + multiple elastic loops for extra storage.

Emergency Response | Backpack Refills



BURN MODULE REFILL
Code: M6001 Case Included

CODE	CONTENTS	QTY
ABD10	AeroBurn Burn Dressing 4in x 4in	2
AF100	AeroForm Conforming Bandage 4in x 4.3yds	1
AGNPF02	AeroGlove Nitrile Gloves Large	2
APD57550BS	AeroPad Low Adherent Dressing 2in x 3in	2
ASC13PUS	AeroSupplies Scissors Blue Lister 5.2in	1
AB0905B	AeroBurn Burn Gel Packet 3.5g	25
ABS6090	AeroBurn PE Burn Sheet 23 1/2 x 35 1/2in	1
ABS2020	AeroBurn PE Burn Sheet 8in x 8in	1
ABS1010	AeroBurn PE Burn Sheet 4in x 4in	1
6001BG	Zippered Bag	1



EYE MODULE REFILL
Code: M6002 Case Included

CODE	CONTENTS	QTY
AEP1S	AeroPad Eye Pads 2in x 3in	4
AF25	AeroForm Conforming Bandage 1in x 4.4yds	1
AF75	AeroForm Conforming Bandage 3in x 4.4yds	1
AGNPF02	AeroGlove Nitrile Gloves Large	2
AEM6PSUS	AeroSupplies Eye Magnet Plastic	1
AW1000US	AeroWash Eyewash .5oz	8
AP312	AeroTape Microporous Paper Tape .5in x 5.5yd	1
6002BG	Eye Wound Zippered Bag	1



GENERAL FIRST AID MODULE
Code: M6006 Case Included

CODE	CONTENTS	QTY
AFP6012S	AeroPlast Premium Fabric X-Wide Strip 1in x 3in	25
AP1012S	AeroPlast Plastic Strip 0.7in x 2.8in	25
AFP312	AeroPlast Premium Fabric Fingertip & Knuckle Dressings	12
AF50	AeroForm Conforming Bandage 2in x 4.3yds	2
AF9	AeroSupplies Forceps Steel 3.5in	1
ASP100	AeroProbe Splinter Probes 1.4in	5
ASC13PUS	AeroSupplies Scissors Blue Lister 5.2in	1
AW8001	AeroWipe Antiseptic Wipe	30
AW7001	Alcohol Swabs 1.25in x 1.25in	10
AGNPF02	AeroGlove Nitrile Gloves Large 1 Pair	1
AT25	AeroPlast Transparent Tape 1in x 5.5yd	1
APD57550BS	AeroPad Low Adherent Dressing 2in x 3in	1
APG1050BSUS	AeroPad Gauze Pad 4in x 4in	2
APG7550BS	AeroPad Gauze Pad 3in x 3in	4
APG550BSUS	AeroPad Gauze Pad 2in x 2in	2
AGG001	AeroGuide First Aid Guide	1
ASTA05US	AeroSoothe Triple Antibiotic 0.5g	25
ACG10US	AeroCleanse Hand Sanitizer 0.9g	20
ACD1223S	AeroWound Combine Dressing 5in x 9in	1
AWF002	AeroWound Finger Dressing	1
M6006_BAG	General First Aid Zippered Bag	1

Emergency Response | Backpack Refills



SPRAINS & FRACTURES MODULE
Code: M6007 Case Included

CODE	CONTENTS	QTY
ATB213	AeroRescue Thermal Blanket 55in x 84in	1
ABN10	AeroBand Triangular Bandage	2
ASP12	AeroSupplies Safety Pins Medium	12
AGNPF02	AeroGlove Nitrile Gloves Large 1 Pair	1
AFH100	AeroForm Conforming Heavyweight Bandage 4in x 4.3yds	1
AFH50	AeroForm Conforming Heavyweight Bandage 2in x 4.3yds	1
All3100US	AeroPlast Instant Ice Pack 4in x 5in	2
AW8001	AeroWipe Antiseptic Wipe	10
AR1060R	AeroSplint Aluminium Foam Splint 4.25in x 24in Rolled	1
M6007_BAG	Sprains & Fractures Zippered Bag	1



STANDARD BLEED CONTROL MODULE
Code: M6008 Case Included

CODE	CONTENTS	QTY
AF100	AeroForm Conforming Bandage 4in x 4.3yds	1
ACD1223S	AeroWound Combine Dressing 5in x 9in	3
AW8001	AeroWipe Antiseptic Wipe 5in x 7in	10
RST103	RapidStop Arterial Tourniquet	1
RSC104	Permanent Marker Black	1
AGNPF02	AeroGlove Nitrile Gloves Large 1 Pair	2
ATB213	AeroRescue Thermal Blanket 55in x 84in	1
AWD100S	AEROWOUND Pressure Dressing 4in x 4in	1
AWD75SS	AEROWOUND Pressure Dressing 3in x 3in	1
AWD50S	AEROWOUND Pressure Dressing 2in x 2in	1
M6008_BAG	Bleed Control Zippered Bag	1



PRO BLEED CONTROL MODULE
Code: M6009 Case Included

CODE	CONTENTS	QTY
AGNPF02	AeroGlove Nitrile Gloves Large 1 Pair	2
ATB213	AeroRescue Thermal Blanket 55in x 84in	1
RST103	RapidStop Arterial Tourniquet	1
RSC104	Permanent Marker Black	1
RSC100	Hemorrhage Control Bandage 4in	1
RSC103	Chest Seal Dressing Vent 2/pack	2
RSC102	Hemostatic Gauze 3in x 4yd	1
USP15	Trauma Shears 5.5in	1
ACD1223S	AeroWound Combine Dressing 5in x 9in	3
AF100	AeroForm Conforming Bandage 4in x 4.3yds	1
AWD100S	AEROWOUND Pressure Dressing 4in x 4in	1
AWD75SS	AEROWOUND Pressure Dressing 3in x 3in	1
AWD50S	AEROWOUND Pressure Dressing 2in x 2in	1
AW8001	AeroWipe Antiseptic Wipe	10
M6008_BAG	Bleed Control Zippered Bag	1

Emergency Response | Backpack Refills



AED / CPR PREP MODULE
Code: M6010 Case Included

CODE	CONTENTS	QTY
AM01BUS	AeroShield CPR Mask & Gloves	1
ARZ001	AeroSupplies Double Sided Surgical Razor	1
AGNPF02	AeroGlove Nitrile Gloves Large 1 Pair	2
AHBB24US	AeroHazard Biohazard Bag 24in x 24in	1
AW8001	AeroWipe Antiseptic Wipe 5in x 7in	1
USP15	AeroSupplies Universal Shears Plastic Handle 5.9in	1
AHPT	AeroHazard Absorbant Towel 10in x 13in	1
M6010_BAG	Yellow Zippered Bag	1



AIRWAY & RESUSCITATION MODULE
Code: M6011 Case Included

CODE	CONTENTS	QTY
36315	Guedel Airway Infant - Pink 40mm	1
36317	Guedel Airway Child - Black 60mm	1
36319	Guedel Airway Small Adult - Green 80mm	1
36320	Guedel Airway Medium Adult - Yellow 90mm	1
36321	Guedel Airway Large Adult - Red 100mm	1
36322	Guedel Airway X-Large Adult - Orange 110mm	1
36275	Nasopharyngeal Airway 20 FR	1
36270	Nasopharyngeal Airway 24 FR	1
36279	Nasopharyngeal Airway 28 FR	3
1250	Lubricating Jelly Packets 2.7g	5
36000	Adult Manual Pulmonary Resuscitator, Mask, Bag, Tubing	1
AGNPF02	AeroGlove Nitrile Gloves Large	2

This sturdy metal wall bracket keeps your backpack visible and ready for emergencies.



Item: Metal Wall Bracket
Code: MWB01



This empty black medical supplies bag is designed to fit securely inside the top compartment of the backpack and can be customized with your own medical tools and supplies.

Item: Modulator Refill - Medication & Supplies (Empty / No Contents)
Code: M6012