

Point-of-Care, Phlebotomy Cart



Features

- Tough, powder-coated, steel exterior, center cabinet
- ClintonClean™ solid HDPE, plastic base and top
- Top has integrated easy-grip handle and 3/4" catch channel to help prevent objects from rolling off the surface and a 4" x 5" integrated cutout to hold items
- ClintonClean™ door with clear PETG panels makes identifying supplies quick and easy
- Easy-clean, white plastic interior sides
- Three 3.5" medical, drawer bins are standard. Choose options for other medical, drawer bin configurations. All medical, drawer bins come with dividers.
- All drawer bins are on rust-free plastic tracks
- Standard 4", hospital style casters, two locking
- 15.75" x 4" glove box holder with storage compartment; three, 15.75" x 4" clear, covered, tilt-lid, side bins; and three, 15.75" x 3" open side bins. All side bins have various size compartments.
- Keyed lock secures* cabinet
- All parts are field replaceable, if ever damaged
- 3-year, limited warranty

NOTE: Supplies pictured not included.

*** Lock not intended to secure narcotics or other restricted items.**

SPECIFICATIONS

Top	25.25" L x 19" W
Base	22.5" L x 19" W"
Height	37.5"
Cabinet	17" L x 14" W
3.5" Blue Bin	15" L x 9" W x 3.5" D
5" Blue Bin	15" L x 9" W x 5" D

Prop 65 Warning



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



MODEL 67-117
Three 3.5" Bins



Open Door



67-117 with option 304
Three 3.5" Bins & One 5" Bin



67-117 with option 305
Four 5" Bins & One 3.5" Bin



67-117 with option 306
Six 3.5" Bins



Top Cutout
4" x 5" Cutout for stabilizing objects. Tops are also available without cutout at no charge. Select "No Cutout Option -611"

Additional Point-of-Care, Phlebotomy Cart Options



-601 Electronic Lock

NOTE: Electronic Lock not intended to secure narcotics or other restricted items.



-602 Five inch Casters
(upgrade from 4 inch)



-611 No 4" x 5" Cutout
—No Charge



-612 Sliding Tray Holder
18.25" W adjusts from 8.25" to 13.5" D



-617 One-Piece, Easy-Clean, Stainless Steel, Push Handle & Small Object Tray

Dividers in all Bins (make up to 6 sections)



DIVIDERS INCLUDED WITH EACH BIN:

- Two Short & One Long

ADDITIONAL DIVIDERS

CHOOSE OPTIONS:

-735 (3.5" DEEP SET)
Three Short & One Long

-750 (5" DEEP SET)
Three Short & One Long